

+1 (732) 230-3003

For inquiries: quote@synbio-tech.com
For assistance: service@synbio-tech.com

SERVICE BROCHURE

A T C G

Oligo
Synthesis

DNA
Library
Synthesis

Gene
Synthesis

Recombinant
Protein
Expression

Antibody
Discovery

SCIENCE
GENES FOR LIFE

QUALITY
GENES FOR LIFE

PEOPLE
GENES FOR LIFE



About Us

Founded in 2013 by a group of passionate synthetic biologists with decades of industry experience and a proven track record of delivering cost-effective DNA solutions, Synbio Technologies is incorporated with a single mission in mind: to empower scientific discoveries and drug innovations by providing the most advanced and cost-effective DNA technology platform. The company takes the lead in establishing an advanced biotechnology transformation and application “GPS” (Genotype-Phenotype-Synotype) platform. The innovative “GPS” Platform follows the biological central dogma, expands the genotype, phenotype, and synotype effectively to achieve the one-stop solution of XNA “design-construction-application”. With a full range of DNA reading (sequencing), DNA writing (synthesis), and DNA editing (engineering) capabilities, Synbio Technologies has delivered satisfactory solutions globally to researchers for various applications, including diagnostic DNA probes, precision medicine, protein production, antibody discovery, vaccine development, novel enzymes, molecular breeding, biofuel implication, etc.

.....▶▶▶▶

Synotype

▶▶▶▶

Design→Build→Test→Learn



Synthetic Biology Dogma



Phenotype

▶▶▶▶

Application

DNA→RNA→Protein



Genotype

◀◀◀◀

Central Dogma

Chapter 1 Syno® C Oligo Synthesis 01

1.1 Standard Oligo Synthesis	03
1.2 Oligo Modification & Labeling	06
1.3 High Quality Oligos for NGS	09
1.4 Oligo Pool Synthesis	10

Chapter 2 Syno® GS Gene Synthesis 11

2.1 Standard Gene Synthesis	13
2.2 Synstrands Gene/DNA Fragments	14
2.3 Pathway Synthesis & Assembly	15
2.4 Small Genome Synthesis & Assembly	16

Chapter 3 DNA Library Synthesis 18

3.1 Precision Mutant Libraries	21
3.2 CRISPR sgRNA Library	25

Chapter 4 Molecular Biology Services 30

4.1 Plasmid DNA Preparation	32
4.2 PCR Cloning & Subcloning	34
4.3 Site-Directed Mutagenesis Service	35

Chapter 5 CRISPR sgRNA Genome Editing Platform 36

5.1 CRISPR sgRNA Design Center	38
5.2 CRISPR System Activity Detection	39
5.3 Ready-to-Use sgRNA Synthesis	41
5.4 CRISPR sgRNA Vector Construction	44
5.5 ssDNA Synthesis	45
5.6 CRISPR Sequencing	46

Chapter 6 Protein Expression and Purification 47

6.1 NG™ Codon Optimization Technology	49
6.2 Guaranteed Bacterial Protein Expression	50
6.3 Yeast Protein Expression System	51
6.4 Baculovirus-Insect Cell Expression System	52
6.5 Mammalian Cell Expression System	53
6.6 Peptide Synthesis	55

Chapter 7 Antibody Services 57

7.1 Monoclonal/Hybridoma Antibody Sequencing	60
7.2 Immune Repertoire Sequencing	62
7.3 <i>De Novo</i> Antibody Design and Production	65
7.4 Antibody Humanization Service	67
7.5 Antibody Affinity Maturation Service	68
7.6 Recombinant Antibody Production	69
7.7 Monoclonal Antibody Preparation	70
7.8 Polyclonal Antibody Preparation	71
7.9 Single Domain Antibody (sdAb) Service	73

Chapter 8 RNA Synthesis 75

8.1 Custom RNA Oligos	77
8.2 RNA Modifications	78
8.3 siRNA Synthesis	79
8.4 miRNA Synthesis	80
8.5 RNA <i>In Vitro</i> Transcription Service	81

Chapter 9 COVID-19 DNA Solutions 83

9.1 Probe Synthesis for COVID-19 Detection	85
9.2 Gene Synthesis & Protein Expression for COVID-19	90

Chapter 10 DNA Data Storage 91

10.1 Your IP, Our DNA	94
-----------------------	----

CHAPTER 1

Syno[®] C Oligo Synthesis

Syno[®] C Oligo Synthesis

Oligo Synthesis

Synbio Technologies's oligo synthesis platform equips state-of-the-art synthesizers and leading professional teams experienced in oligo synthesis and modification technologies. This allows us to provide high-quality oligo products to scientific researchers and industrial customers around the world.


By leveraging our advanced and innovative manufacturing process, Synbio Technologies offers a wide range of oligo synthesis options including standard oligos, long oligos, ultra-high purity oligos, degenerate primers, aptamers, fluorescent probes, oligo pool synthesis, etc. In addition, we also provide high-throughput oligo synthesis to customers with large-scale oligo projects.

Competitive Advantages



- ▶ Certified with the highest manufacturing standards (ISO9001 & ISO13485), large-scale sequencing that verifies the average mutation rate to be lower than 1/1000
- ▶ Various highly customizable options for purification, synthesis specifications, and delivery form.
- ▶ Qualified and available for multiple applications such as molecular diagnostics, NGS, functional genomics, genome engineering, drug development, protein/antibody function, etc.

Oligo Synthesis Services



1.1 Standard Oligo Synthesis

Synbio Technologies is always competent to offer tailor-made approaches to accomplish every oligo synthesis request.

Competitive Advantages

High Quality

We can achieve extremely low mutation and error rates in accordance with stringent oligo synthesis quality control standards

Highly-Customizable

Flexible synthesis scales are available. Four alternative purification options include: DSL, HPLC, PAGE Plus, and PAGE

Technical Support

Professional teams experienced in oligo synthesis and various modifications offer support throughout the production progress

Cost Effective

Competitive prices and affordable services

Purification Methods

Purification Method	Characteristics	Purity	Applications
DSL	Effectively removes salt and meets most experimental needs	>75%	It is used for PCR cloning, site-directed mutagenesis, gene synthesis, Sanger sequencing, etc.
PAGE-Plus	Purification effect is better than DSL, especially suited for oligos longer than 59 nt and meets most experimental needs	>85%	Commonly used in sequencing, multiplex & quantitative PCR cloning, site-directed mutagenesis, gene synthesis, etc.
HPLC	Purification effect has high degree of automation and is good for modified oligos and effectively removing small fragments	>90%	Often used in modified oligos, probes, cloning, etc.
PAGE	Purification effect is good for removing small fragments and is generally more suitable for oligos of more than 59 nt or for high quality requirements	>95%	Commonly used in sequencing, gene synthesis, site-directed mutagenesis, PCR cloning, etc.

◆ Standard Oligos

Synbio Technologies provides fast and convenient standard DNA oligo synthesis services with no design fee that are suitable for more than 95% of PCR projects. Large-scale and bulk custom DNA oligos can be ordered with the same high quality standards, competitive prices, and fast turnaround.

	Purification Method	Synthesis Scale	10 nm	25 nm	50 nmol	100 nmol	200 nmol	1 μmol
11-59nt	DSL	Final Yield	1 OD	2 OD	5 OD	10 OD	15 OD	50 OD
	PAGE Plus	Final Yield	–	–	2 OD	5 OD	10 OD	30 OD
	PAGE	Final Yield	–	–	–	2 OD	5 OD	10 OD
	HPLC	Final Yield	–	–	–	2 OD	5 OD	10 OD
60-90nt	DSL	Final Yield	1 OD	2 OD	5 OD	10 OD	15 OD	50 OD
	PAGE Plus	Final Yield	–	–	2 OD	5 OD	10 OD	30 OD
	PAGE	Final Yield	–	–	–	2 OD	5 OD	10 OD
	HPLC	Final Yield	–	–	–	2 OD	5 OD	10 OD

* Longer oligos, 96/384 plate batch oligos, or other large scale synthesis are available, please contact us at quote@synbio-tech.com.

* Please refer to the purification table for the final output corresponding to each synthesis scale.

◆ Degenerate Primers

Genetic codes have the tendency to degenerate, which lead to the use of degenerate primers in PCR experiments to exhibit specific genes. Synbio Technologies provides both the design and synthesis of these degenerate primers, as well as custom degenerate oligo pool synthesis upon request.

Symbol	Base Type	Standard Ratio
R	A, G	50% : 50%
Y	C, T	50% : 50%
M	A, C	50% : 50%
K	G, T	50% : 50%
S	C, G	50% : 50%
W	A, T	50% : 50%
H	A, C, T	33% : 33% : 33%
B	C, G, T	33% : 33% : 33%
V	A, C, G	33% : 33% : 33%
D	A, G, T	33% : 33% : 33%
N	A, C, G, T	25% : 25% : 25% : 25%

* Non-standard base proportion should be declared clearly, for example, N (20%A: 30%C: 40%G: 10%T).

CHAPTER 1

Syno[®]C Oligo Synthesis



◆ Trimer Oligos

Utilized for certain projects sensitive to codon bias. Synbio Technologies can significantly reduce the deviation inherent in the synthesis procedures and realize accurate proportional control. The basic principle is to introduce 20 amino acids into the specific sites of proteins by the chemically synthesized trimer phosphoramidites, reduce codon redundancy, avoid codon shift mutations, and stop codon insertion.

Trimer Oligo Length	The Number of Trimer Sites	Purification Methods	Turnaround Time	Deliverables
≤120 nt	≤20	PAGE/HPLC	2-6 weeks	<ul style="list-style-type: none">• Trimer oligos• COA• NGS analysis report (Optional)

◆ Long Oligos

Synbio Technologies’s oligo synthesis platform utilizes an optimized ultra-long chain design and synthesis program to meticulously monitor and control the synthesis process. This process allows us to generate long, high-quality oligos (90-129 nt) and even ultra-long oligos (130-150 nt) according to individual requests.

Length (nt)	Quantity	Purification Methods	Turnaround Time (Business Days)
90-109	1 OD	PAGE / Dual-PAGE	5-6
110-129	1 OD	PAGE / Dual-PAGE	5-7
130-150	0.5 OD	PAGE / Dual-PAGE	5-8

Standard Deliverables

1. One tube lyophilized DNA (transparent or dark) or 96-well and 384-well plate
2. Certificate of analysis (COA) including sequencing information, OD, Tm, etc.

1.2 Oligo Modification & Labeling

Synbio Technologies provides comprehensive modification and labeling services to meet our customers’ specified research requirements. Our stringent quality control guidelines ensure the synthesis of high-quality modified oligos and that they are subsequently delivered in a timely manner. We can satisfy various research needs in biology, diagnostics, personal medicine, drug discovery, etc.

Competitive Advantages

High Synthesis Efficiency

Large-scale sequencing verifies that the average mutation rate will be lower than 1/1000

Strict Quality Control

Stringent ISO 9001- & ISO13485-certified processing for oligo synthesis, meticulously monitored with advanced traceability

Comprehensive Modifications and Labeling

Professional teams experienced in oligo synthesis and various modifications offer support throughout the production progress

Service Specifications

◆ General Modifications

We provide a series of general modifications with high purity and competitive prices. We accomplish this by focusing on modification requests with special chemical classifications.

Classifications	Modifications	Classifications	Modifications	Purifications
Modified Bases	dU	AMCA	5' AMCA	HPLC/PAGE
	dI		3' AMCA	
	5-Methyl dC		dT-AMCA	
DBCO	5' DBCO	Biotin	5'Biotin	
	3' DBCO		3'Biotin	
	dT-DBCO		dT-Biotin	
Amino Linkers	5'Amino C6	Fluoro Bases	2'-F-2'-dA	
	dT-Amino C6		2'-F-2'-dC	
	3'Amino C7		2'-F-2'-dG	
	5'Amino C12		2'-F-2'-dU	
Digoxigenin	5'Digoxigenin	Spacer	dSpacer	
	3'Digoxigenin		Spacer C3	
	dT-Digoxigenin		Spacer C9	
Locked Nucleic Acid	LNA (A/G/C/U)		Spacer C18	

Classification	Modifications	Classification	Modifications	Purification
Phosphorylation	5' PO4	Hydrosulphonyl	5' HS-C6	HPLC/PAGE
	3' PO4		3' HS-C3	
Phosphorothioates	SPO3	—	—	

* If you need other purification methods, please contact us at quote@synbio-tech.com.

◆ Single Labeling

Synbio Technologies provides multiple fluorescent groups at 3', 5', or other designated sites to comply with our customers' requests. These include FAM, HEX, ROX, BHQ, Texas Red, etc.

Dye/Quenchers	3' -	5' -	Designated Sites	Purification
FAM	√	√	√	HPLC/PAGE
TET	√	√	√	
JOE	√	√	√	
HEX	√	√	√	
SIMA	√	√	√	
TAMRA	√	√	√	
ROX	√	√	√	
Cy5	√	√	√	
Cy3	√	√	√	
Quasar670	√	√		
Quasar570	√	√		
DyLight547		√		
DyLight647		√		
Texas Red		√		
BHQ-1	√			
BHQ-2	√			
BHQ-3	√			
Eclipse	√			
DABCYL	√			

* If you need other purification methods, please contact us at quote@synbio-tech.com.

◆ Dual/Multiple Labeling

Aside from single fluorescent groups, Synbio Technologies also offers dual/multiple labeling services including TAMRA, BHQ, DABCYL, ELCISPE, etc. The most common dyes and quenchers are listed below:

Dyes	Quenchers	Dyes	Quenchers	Purification
5'6-FAM	3'TAMRA	5'HEX	3'TAMRA	HPLC/PAGE
	3'BHQ-1		3'BHQ-1	
	3'Eclipse		3'BHQ-2	
	3'DABCYL		3'DABCYL	
5'TAMRA	3'Eclipse	5'ROX	3'Eclipse	
	3'BHQ-1		3'BHQ-1	
	3'BHQ-2		3'BHQ-2	
	3'DABCYL		3'DABCYL	
5'CY5	3'BHQ-1	5'CY3	3'BHQ-2	
	3'BHQ-2	5'Texas Red	3'BHQ-2	
	3'BHQ-3	5'JOE	3'Eclipse	
	3'DABCYL		3'BHQ-1	
5'TET	3'TAMRA		3'BHQ-2	
	3'BHQ-2		3'DABCYL	
	3'DABCYL		3'TAMRA	

* If you need other purification methods, please contact us at quote@synbio-tech.com.

Standard Deliverables

1. One tube lyophilized DNA (transparent or dark) or 96-well and 384-well plate
2. Certificate of analysis (COA) including sequencing information, OD, Tm, etc.

1.3 High Quality Oligos for NGS

Next Generation Sequencing (NGS) is a high-throughput technology that enables the rapid sequencing of millions of DNA strands in parallel. With its ultra-high throughput capabilities, scalability, and speed, NGS has revolutionized genomics and biological sciences. With NGS, researchers can perform a variety of applications as well as study biological systems beyond the capacity of traditional sequencing technologies.

For most of the NGS platforms available on the market, there is a large demand for high quality oligos for successful sequencing. The quality of the oligos are crucial because any contamination can lead to strong background or splicing errors. These errors may seriously affect the overall sequencing results while also wasting valuable time and effort in the process. Synbio Technologies's Syno[®] C oligo synthesis platform allows us to avoid these potential errors while delivering custom, high-quality oligos for NGS applications.

Competitive Advantages

01

Industry-Leading Coupling Efficiency
Reduce the dimer generation of synthetic oligos

02

Exclusive Production Processes
Standardized and stringent quality control ensures batch-to-batch consistency

03

Low Levels of the Cross-Contamination
Strict purification standards generate higher quality products, often with purity ≥90%

04

Flexible Oligo Specifications
Customizable formulation, mixing options, and documentation are available

Custom High Quality Oligos for NGS Services

Our NGS-grade oligos are recommended for a variety of requests independent of the instrument or technology. Customized synthesis, modification, purification, and delivery are provided to meet all our customers' needs.

- Syno[®] Adapters
 - Syno[®] Hybrid Capture Probes
- Syno[®] Blocking Oligos
 - Syno[®] UMI
- Amplification Oligos
 - Sequencing Oligos

Specifications	Length (bp)	Turnaround Time (Business Days)	Deliverables
Indexed Adapters	60-70	5-7	Products in tube/plate/chip, Comprehensive QC report (HPLC, MS)
Capture Probes (DNA oligo – 120 nt biotin, Oligo pools, Multiplex PCR oligos)	10-170	Inquire	
Blocking Oligos	~60	5-7	
Other NGS oligos	10-150	5-7	

1.4 Oligo Pool Synthesis

Many aspects of biological research require a large amount of primers, such as genetic library construction and target enrichment/capture in NGS. Conventional approaches to synthesizing hundreds of thousands of primers at a time are ineffective and time consuming. To circumvent this issue, Synbio Technologies provides a solution to high-throughput DNA synthesis with our industry-leading technology.

Synbio Technologies's Syno[®] HT synthesis platform can synthesize 890,000 oligos simultaneously on a semiconductor chip with single oligo lengths up to 300 nt. High synthesis efficiency can ensure high uniformity of the oligo pool, improving the efficiency and success rate of downstream high-throughput screening while perfectly matching customers' experimental needs and applications.

Competitive Advantages

High Throughput

Up to 890,000 sequences are synthesized in parallel at a single time

High Accuracy

Error rate as low as 1:2,000/ nt

High Coverage

100% coverage of the target sequences

High Homogeneity

More than 90% of oligonucleotides vary within 2.5 times

Cost-effective

Multiple streamlined services available at minimal cost

Service Specifications

Length (nt)	Turnaround Time (Business Days)
20-120	10-15
121-150	10-15
151-200	10-15
201-250	15-20
251-300	20-25

Multiple Downstream Applications

Synbio Technologies's oligo pool services offer our customers industry-leading high-throughput synthesis techniques as well as high-quality constructs. With the synthesized oligo pools, our customers can then streamline the various downstream applications within their research.

- Generate native DNA sequences
- Create *de novo* DNA sequences
- Build genes, chassis, operons, pathways, and genomes
- Build DNA variant libraries
- Improve the features of a protein
- Test all orthologs of a gene

- Optimize any antibody through affinity maturation
- sgRNA screening library
- shRNA library
- NGS hybridization capture probe library
- FISH probe

CHAPTER 2

Syno[®] GS Gene Synthesis

Gene Synthesis

While Synbio Technologies offers many different services, our main focus is on DNA synthesis. With our professional research team and comprehensive Syno[®] systems, we can provide global customers with various types of DNA synthesis services, such as gene synthesis, gene fragment synthesis, pathway synthesis, small-genome synthesis, etc. We routinely synthesize over 3 million base pairs of DNA sequences each month for a wide range of applications including protein production, antibody discovery, genetically engineered vaccine, molecular breeding, biofuel, and many more.



- ▶ 100% sequence accuracy guaranteed with success rate >99.5%
- ▶ AI-aided sequence analyses and proprietary NG[™] codon optimization
- ▶ Industry's most competitive rate
- ▶ FREE to subclone into most commercial vectors
- ▶ FREE to deposit your vector for future use

Gene Synthesis Services




- Standard Gene Synthesis
- High Throughput Gene Synthesis
- NG[™] Codon Optimization
- Vector Construction
- Small Genome Synthesis & Assembly
- Pathway Synthesis & Assembly
- Synstrands Gene Fragment Synthesis

Syno[®] GS
Gene Synthesis

2.1 Standard Gene Synthesis

According to the DNA or amino acid sequences that customers provide us, our scientists will design the oligos and the procedure of fragment assembly and cloning within 30 minutes. By using our automated system that guarantees efficient and error-free synthesis, our customers will be provided with the highest quality output. Sanger sequencing will be conducted to verify 100% sequence accuracy.




◆ Service Specifications

The following turnaround times only apply to non-complex sequences. Gene products will be cloned into our default pUC57-Amp or pUC57-Kan vectors with multiple cloning sites. We can also clone genes into customers' supplied vectors upon request.

Gene Length	Turnaround Time (Business Days)	Standard Deliverables	Optional Services
< 250 bp	5-10	• 2-5 µg lyophilized plasmid DNA • Sequencing chromatogram • Certificate of analysis	• Free codon optimization • Free vector deposit
250-1,500 bp	5-10		
1,501-3,000 bp	10-15		
3,001-4,500 bp	15-20		
4,501-6,000 bp	20-25		
> 6,000 bp	Quote		

* We offer discounts for bulk gene synthesis service. Please contact quote@synbio-tech.com for additional information.


Synbio Technologies independently developed a CI sequence difficulty index system, which can facilitate the accurate, fast, and high-quality synthesis of any sequence.




2.2 Synstrands Gene/DNA Fragments

Synbio Technologies supplies our customers with double-stranded, sequence-verified Synstrands gene fragments in only a few business days. These fragments offer affordable and accessible construction or modification of genetic sequences. The various applications include antibody research and CRISPR-mediated gene/genome editing, being used as qPCR standards, etc. Synstrands gene fragments are available in two delivery forms: eppendorf tubes or plates.

Competitive Advantages



Process of the Service



Service Specifications

◆ Normal Sequences

Length (bp)	Turnaround Time (Business Days)
<300	8-10
300-1,800	8-10
1,801-3,200	10-12

2.3 Pathway Synthesis & Assembly


The synthesis and construction of metabolic pathways, as well as diverse DNA assembly technologies, play fundamental roles in microbiology, biochemistry, and many other relevant fields. With their rapid development, the price of large DNA fragment or plasmid synthesis is decreasing year by year. Based on our Syno[®] gene synthesis tools and excellent long DNA assembly capability, Synbio Technologies is confident in providing professional technical support for metabolic pathway synthesis and library assembly.

Services Specifications

Services	Details	Turnaround Time	Delivered Products
Metabolic Pathway Synthesis	Customized Metabolic Pathway Synthesis Services	<10 Kb, 4-8 weeks >10 Kb, Inquiry	Plasmids with synthesized metabolic pathway and glycerol bacteria
Metabolic Pathway Library Assembly	Customized Metabolic Pathway Library Assembly Services	Varies depending on the pathway requested	Plasmid library after assembly and glycerol bacterial library, or monoclonal plasmids and glycerol bacteria

Case Study

The promoter library's diversity determines the final diversity of the constructed library. The theoretical library capacity is 100, with 100 times coverage, leading to the production of 10,000 clones. Our PCR test of the bacterial solution showed 78% positive and 100% accuracy.



2.4 Small Genome Synthesis & Assembly

The rapid development of synthetic biology has allowed for researchers to place greater emphasis on powerful concepts and techniques within the field itself. These techniques include the comprehensive analysis of complex pathways and the re-design of natural biological systems, such as Gibson Assembly. Currently, Gibson Assembly is the most popular mechanism for achieving directed cloning and *in vitro* multiple-segment assembly. However, Gibson Assembly is limited to plasmids around 10 Kb or less in length, which largely restricts its applications within certain aspects of synthetic biology. Several researchers discovered that assembly in yeast has the potential to conveniently assemble large DNA segments with ease. This is because yeast can connect multiple DNA segments automatically without DNA polymerases and ligases. This important aspect of yeast made it possible to circumvent the length restriction associated with Gibson Assembly, and efficiently construct a plasmid over 10 kb in length. For this reason, yeast homologous recombination technologies are greatly beneficial for efficient and effective assembly of long DNA fragments and genomes.

Competitive Advantages

01

Powerful gene synthesis platforms provide one-stop solutions for yeast assembly platforms such as gene synthesis, assembly, and sequencing

02

Synthesis and assembly of long DNA fragments with high fidelity and speed. Synbio Technologies can deliver single gene sequences up to 150 Kb, guaranteeing reliable assembly of long DNA fragments


Services Specifications

Services	Details	Turnaround Time	Deliverables
Synthesis of gene clusters and small genomes	Customized service for metabolic pathway synthesis	<10 Kb, 4-8 weeks >10 Kb, Inquiry	Synthethized gene clusters, plasmids withsmall genomes, and glycerol bacteria

CHAPTER 2

Syno[®] GS Gene Synthesis

By utilizing yeast homologous recombination technology, Synbio Technologies is able to provide one-stop services for long gene synthesis/assembly and gene cluster/small genome synthesis and assembly.




DNA Library
Synthesis

CHAPTER 3

DNA Library Synthesis

DNA Library Synthesis




Synbio Technologies offers high quality DNA libraries of any size and complexity with an efficient time frame and at the industry's most competitive rate. Our synthetic DNA libraries are designed by an experienced bioinformatics team to facilitate various research applications, including functional genetics, directed evolution, protein engineering, antibody studies, etc.



- ▶ Proprietary synthesis platforms guarantee the successfully synthesis of multiple libraries
- ▶ High-quality libraries can be accurately synthesized with capacity from 10⁶ to 10⁹, coverage rate>99%, uniformity index<10, NGS verification accuracy>95%
- ▶ Flexible deliverables: oligo pools, PCR fragments, pooled plasmids, plasmid arrays, etc

Library Types

Library Types	Descriptions	Schema
Alanine Scan Library	Research on amino acid residues critical to protein function, interaction, and shape by replacing an individual amino acid with alanine at every position.	<div><div>F</div><div>A</div><div>D</div><div>F</div><div>E</div><div>K</div><div>E</div><div>D</div><div>F</div><div>I</div></div> <div><div>F</div><div>I</div><div>A</div><div>F</div><div>E</div><div>K</div><div>E</div><div>D</div><div>F</div><div>I</div></div> <div><div>F</div><div>I</div><div>D</div><div>A</div><div>E</div><div>K</div><div>E</div><div>D</div><div>F</div><div>I</div></div> <div><div>F</div><div>I</div><div>D</div><div>F</div><div>A</div><div>K</div><div>E</div><div>D</div><div>F</div><div>I</div></div> <div><div>F</div><div>I</div><div>D</div><div>F</div><div>E</div><div>A</div><div>E</div><div>D</div><div>F</div><div>I</div></div>

Library Types	Descriptions	Schema
Antibody Library	Create a high diversity synthetic antibody library or improve existing antibody functionality such as specificity, immunogenicity, affinity, expression, and/or aggregation.	
Modular Substitution Library	Build new genetic modules, regions, or pathways from DNA parts including regulatory elements and/or coding sequences.	
CRISPR Library	Customized CRISPR library and human genome CRISPR knockout library meet the screening needs of different customers.	

3.1 Precision Mutant Libraries

In vitro synthesis of precision mutant libraries are one of the common methods in molecular directed evolution. Researchers can change specific protein sites or DNA regulatory regions according to different purposes, while controlling the diversity of sequences and improving screening efficiency. With more than ten years of experience in precision mutant library synthesis, Synbio Technologies provides customers around the world with efficient and accurate precision mutant library synthesis services. We are committed to providing customers with efficient and cost-effective solutions for downstream research, such as protein characteristic improvement, protein directed evolution, antibody sequence optimization, affinity maturation, target discovery, etc.

Competitive Advantages

Accurate
Sequence
Control

High
Synthesis
Efficiency

Cost-Effective
Solution

Highly
Customized

Services Specifications

Library Type	Complexity Analysis	Diversity Analysis by Sequencing	Turnaround Time (Business Days)	Deliverables
Site-Saturation Mutagenesis Library, eg.NNK Library	√	√	20-40	• 50-1000 µg plasmids • COA files
Combinatorial Mutagenesis Library, eg Trimer Library	√	√	20-40	

Degenerate Mutation Library

Genetic codons have degeneracy, that is, an amino acid can be encoded by more than one triplet codon. Using the degeneracy of genetic code can not only maintain the diversity of amino acid sequences, but also reduce the redundancy of codons and improve screening efficiency.

Competitive Advantages

01

Quick & efficient product delivery to promote downstream screening

02

The synthesis technology is stable and the quality control is strict

03

100% coverage of 20 amino acids fully ensures the diversity of the library

Case Study

The NNK library is a classic degenerate mutation library. The NNK degenerate primers used contain 32 (4x4x2) codon combinations (N= A/C/G/T, K= G/T), covering all 20 amino acids, and can saturate amino acids at any site.




Fig.1 The proportion of each codon contained in the five NNK sites in the NNK library

CHAPTER 3

DNA Library Synthesis



Trimer library

Primers synthesized with a single nucleotide will inevitably produce redundant mutations and termination codons at the mutation site. The trimer phosphoramidites formed by connecting the three nucleotides in a fixed order can correspond to a single amino acid. It can not only introduce the mixture of all 20 amino acids at any position of the mutated sequence, but also avoids the problems of codon offset and the incorporation of termination codons.

Precise mutant libraries designed with the trimer mixture as the element, also known as trimer libraries, can ensure that the amino acids at the mutation site are composed according to the designed ratio to meet the amino acid diversity and codon preference of the specific site. Trimer libraries are an efficient solution for complex libraries and can further improve the accuracy of library screening.

Competitive Advantages

- 01

Cost-effective solutions for complex libraries
- 02

Reduces difficulty of downstream screening
- 03

High synthesis efficiency and accurate sequence control

Case Study 1

Trimer libraries can be used for site saturation mutation to distribute 20 amino acids in equal proportion, further reducing codon redundancy and removing termination codons in the library.




Fig.1 Trimer libraries are also a saturated mutation library with equal proportion distribution of each amino acid. The proportion of each codon contained in the five sites.

Case Study 2

Trimer libraries can also customize each site to be composed of different ratios of amino acids. Among the four Trimer sites, for example, the X1 Trimer site is A (80%) : G (3%) : R (3%) : S (3%) : T (8%) : V (3%)




Fig.2 Types and proportions of custom amino acids




Fig.3 Ratio between actual and theoretical proportion of each amino acid

3.2 CRISPR sgRNA Library

CRISPR-Cas9 technology has become a powerful tool for forward genetic screening. Screening libraries of various functions can be constructed by CRISPR, and the genes related to this function are then identified through functional screening and enrichment, PCR amplification, and deep sequencing analysis.

Synbio Technologies can provide one-stop services such as sgRNA design, oligo pool synthesis, sgRNA library construction, NGS validation, virus packaging, and bioinformatics analysis. With our professional CRISPR-Cas9 sgRNA library system, we can construct sgRNA libraries in a fast and efficient manner.

Competitive Advantages

01

Advanced sgRNA Design
Design of sgRNA of dozens of species

02

High Efficiency
Designed sgRNA has a low off-target rate and shows strong editing efficiency


03

Proprietary DNA Synthesis Platform
Fast and efficient sgRNA library construction

04

One-Stop Solutions
Customized services ranging from sgRNA design, oligo pool synthesis, sgRNA library construction to lentivirus packaging and NGS sequencing

CRISPR Cas9 sgRNA Library Process



Service Specifications

Service Steps	Services Specifications	Turnaround Time
sgRNA Design	sgRNA sequences were designed for specific genes. Each gene will be designed with 3-6 sgRNAs.	1-2 weeks
Oligo Pool Synthesis	All designed sgRNAs and negative controls will be synthesized on one chip	3-4 weeks
sgRNA Library Construction	Cloning these sgRNAs on the target vector	2-3 weeks
Plasmid Purification	Plasmid preparation to meet the shipping need	1 week
Library Quality Evaluation	1. Library coverage is greater than 100 times 2. Cloning accuracy rate is more than 80% 3. Library accuracy rate of NGS verification is greater than 90% 4. Library uniformity index is less than 10	2-3 weeks
Total		9-13 weeks

Syno® Human Genome CRISPR Knockout Libraries


Synbio Technologies focuses on the research of CRISPR gene editing technology, continuously working to increase the number of services and products offered. Our Syno® Human Genome CRISPR Knockout Libraries select sgRNA sequences with a high score according to the phenotypes in the previously published CRISPR screening recorded in the Genome CRISPR database, based on the criteria of high target activity and low off-target efficiency combined with a variety of different design principles. According to recently published scoring rules, Syno® CRISPR libraries scored better than GeCKOv2 libraries and randomly selected published sgRNA samples.

Syno® Human Genome CRISPR Knockout Libraries consist of three independent sub libraries: A, B, and C. Each gene in Library A and B contains four sgRNA sequences for 18,913 and 18,334 protein coding genes respectively. Sub library A contains the best ranked sgRNA according to the design criteria to achieve high-quality screening under low library coverage. Sub Library B contains a second layer of sgRNA, which can be used to supplement sub library A when higher sgRNA coverage is required.


By integrating the standardized resources of large-scale CRISPR screening and multiple efficiency indicators, sub Library C evaluates and sorts more than 300,000 unique sgRNAs in the CRISPR library, and designs an optimized human minimum genome knockout library. For 18,761 genes, 2 optimized sgRNAs are designed for each gene, in total, there were 37,522 genes targeted and 200 non-targeted sgRNAs. The size of sub Library C is smaller than the currently published human genome CRISPR library. Although the sub Library C is small, it does not affect the identification of necessary genes. It is mostly used for CRISPR screening of cancer cell models under complex conditions.

Competitive Advantages

1. Accurate sgRNA Design: sgRNA sequences with high target activity and low off-target rate are designed by using effective sgRNA design models and deep learning algorithms.



2. Optimized Vector: LentiCRISPRv2 vector with Puro drug screening function can be used for lentivirus packaging and improve transfection efficiency.



3. Cost Saving: High editing efficiency CRISPR library reduces the difficulty of later screening and further reduces research cost.

4. Strict Quality Control: The coverage of the constructed library is > 99.9%, and the uniformity is < 10.

CHAPTER 3

DNA Library Synthesis

Service Specifications

Catalog	Product Name	Sequence Information	Vector Information	Specifications
STCRI211101	Syno® Human Genome CRISPR Knockout Libraries	Sub library A 18,913 genes are targeted, 74,552 sgRNA sequences and 435 control (non-targeted) sgRNA sequences are designed, with a total of 74,987 sgRNA sequences.	lentiCRISPR v2 (with Puro mark)	200 µg
STCRI211102		Sub Library B 18,334 genes are targeted, 71,048 sgRNA sequences and 435 control (non-targeted) sgRNA sequences are designed, with a total of 71,483 sgRNA sequences.	lentiCRISPR v2 (with Puro mark)	200 µg
STCRI211103		Sub Library C 18,761 genes are targeted, and 2 sgRNA sequences are designed for each gene, a total of 37,522 sgRNA sequences and 200 control (non targeted) sgRNA sequences.	lentiCRISPR v2 (with Puro mark)	200 µg

Molecular Biology
Services


CHAPTER 4

Molecular Biology Services



Molecular biology is the cornerstone of many life science research fields, such as biopharmaceutical and human biology. For years, molecular biology technologies have been applied to various fields of biological research and development including molecular mechanisms, drug interactions, targets and effects, etc.

Synbio Technologies provides a wide range of molecular biology services to meet our customers' specific requests and research needs. By relying on revolutionary breakthroughs in DNA synthesis and engineering technologies, our comprehensive services range from plasmid DNA preparation and site-directed mutagenesis to vector construction. These services have allowed Synbio Technologies to become one of the premier companies within the molecular biology technology industry. Compared to the "do it yourself" process within molecular biology services, Synbio Technologies is interested in saving our customers from consuming too much time and energy in laborious operations with fluctuating, sometimes unreliable, results. Our molecular biology services provide our customers with cost-effective and efficient final products. Synbio Technologies can decrease your overall research time and save your lab on both cost and headaches. This process allows our customers to focus on various other research topics in your field of interest, while Synbio Technologies prepares your orders.



Plasmid DNA Preparation

- Fast delivery: starting at one day
- Low endotoxin level: <5 EU/mg
- Diversity of scale: microgram to gram
- ISO 9001 quality management



PCR Cloning & Subcloning

- 100% sequence accuracy
- Customizable service: no limits
- Save your time & budget
- Professional technical team



Site-Directed Mutagenesis


- Short TAT: starting at 5 business days
- Vast region: ≤40 bp region as one mutation
- Precise mutation: error-free
- Large-scale production, such as performing alanine scan library

4.1 Plasmid DNA Preparation


For years, Synbio Technologies has provided comprehensive plasmid DNA preparation services to our customers from research institutes in both biotechnology and biopharmaceutical industries.

Synbio Technologies conducts comprehensive quality control of our plasmid products to provide you with aseptic plasmids with no RNA pollution or genome pollution. All of our manufactured plasmids are free of animal-derived materials and can contain low levels of endotoxin (<100EU/mg, <30EU/mg, <5EU/mg, on request). According to our customers' different application needs, standard plasmid DNA preparation services are divided into two levels, the research level and the transfection level. The plasmid DNA preparation process is designed to meet our customers' diverse downstream applications such as transfection, antibody preparation, vaccines, gene therapy research, etc.

Competitive Advantages



Experimental Procedure



Service Specifications

	Research Grade	Transfection Grade	TAT	Deliverables
Midi (100 µg)	Endotoxin: NA	Endotoxin: < 100EU/mg Endotoxin: < 30EU/mg Endotoxin: < 5EU/mg	2-4 BDs	• Prepared plasmid DNA • Certificate of analysis (COA)
Maxi (200-500 µg)			2-4 BDs	
Mega (1 -10 mg)			3-5 BDs	
Giga (>10 mg)			Starting from 5 BDs	

This form is used for high-copy plasmids. Fore low-copy plasmids, please inquire quote@synbio-tech.com.

Please Note:

1. Sanger sequencing and restriction enzyme analysis are available per request.
2. If you wish to provide the plasmid template, you will need to prepare a sample over 1µg (diluted with ddH₂O or TE buffer), colonies (fresh), or bacteria stored in glycerol.
3. Synbio Technologies also provides upstream services of plasmid preparation such as gene synthesis and subcloning into both commercial and custom vectors. This allows us to offer a unique approach to satisfy all our customers' requests.

Quality Control

Synbio Technologies’s plasmid preparation service focuses on quality. We have strict quality control which can increase quality, stability, and reduce batch difference to satisfy customers’ high-quality needs.

Quality Control	Standard Range	Measurement Technique
Appearance	Clear, transparent	Visual inspection
A260/280 Analysis	1.80-2.00	UV spectrophotometer
Concentration	Up to 5mg/ml, default 1±0.05 mg/ml	UV spectrophotometer
Supercoil Percentage Analysis	>95%	Density scan on agarose gel
RNA Residue Analysis	Not detectable	Visual inspection on agarose gel
Genomic DNA Analysis	Not detectable	Visual inspection on agarose gel
Sterility Assurance	No clone growth on LB plate >48 h	Incubate 0.1 mg of product on LB plate for >48 h
Restriction Enzyme Analysis (optional)	As expected, no other minor band	Incubate with enzyme and analyze on agarose gel
Sequence Verification (optional)	Correct	DNA sequencing
Endotoxin Analysis (optional)	<100EU/mg, <30EU/mg, <5EU/mg, on request	Endotoxin analysis assay

4.2 PCR Cloning & Subcloning

PCR cloning and subcloning are two main approaches to amplifying DNA sequences. Synbio Technologies is confident in our ability to clone any target gene into any requested site of any vector on the basis of our Syno® Platform. The combination of the DNA synthesis tool and the cloning technology allows us to insert the requested genetic sequence into any customer-specific vector for cloning.

During PCR cloning, Synbio Technologies offers sequence verification before introducing target gene materials into selected vectors. This sequence verification allows us to verify the sequence to be 100% accurate to the customer’s request before amplification. According to the model clone and sequence from customers, we will design the necessary primers, clone the PCR product of interest into a new vector, and finally verify it by Sanger sequencing. This allows us to tailor-make a system to supply our customers with the exact specifications necessary to fit their particular research interest.

Competitive Advantages

100%
Sequence
Accuracy

We will quickly and accurately synthesize your requested sequence. The generated sequence will then be verified by Sanger sequencing to guarantee 100% sequence accuracy

Customizable
Subcloning

Guarantee successful cloning at any site in any vector system of interest, both commercial and custom, with special amplified primer design

Professional
Technical
Team

Synbio Technologies's professional technical support offers you proactive communications on your project's progress

Service Specifications

Service Type	Fragment Length	Turnaround Time (Business Days)	Deliverables
PCR Cloning & Subcloning	<1 Kb	5-10	• 2-5 µg Lyophilized plasmid DNA • Sequencing chromatogram • Certificate of analysis (COA)
	1-2 Kb	5-10	
	2-3 Kb	5-10	
	>3 Kb	Inquire	

* Note: The turnaround time applies to non-complex sequences. Please contact quote@synbio-tech.com to quote for the complex sequences.

CHAPTER 4

Molecular Biology Services

4.3 Site-Directed Mutagenesis Service

If the specific site on your gene or vector is not what you want or you need to carry out functional studies of genes or elements, you may need to do a point mutation on a specific site or multiple sites. Site-directed Mutagenesis is a very useful tool for gene research, which can efficiently change the character and characterization of the target protein. Generally, we introduce site-directed mutagenesis through PCR, which can include insertion, deletion, or point mutation. Site-directed Mutagenesis technology has been widely used when studying protein functional site structure, enzyme activity optimization, DNA component function or component interaction, and gene therapy.

Competitive Advantages

Short Turnaround Time

Common sequences can be delivered to you within 5 business days

Vast Operating Region


All mutations found within a 40 bp region will be considered as one whole mutation; This helps our customers save on the overall cost of their project

One-stop Solutions

We provide one-stop solution services, including: gene synthesis vector construction and protein expression and purification, which will certainly improve efficiency

Service Specifications

Cloned Fragment Length	Point Mutation Amounts	TAT (Business Days)	Deliverables
<2 Kb	1	10	• 2-5 µg lyophilized plasmid DNA • Sequencing chromatogram • Certificate of analysis (COA)
	2	10	
	3	10	
2-3 Kb	1	13	
	2	13	
	3	13	
3-5 Kb	1	16	
	2	16	
	3	16	



CRISPR sgRNA

Genome Editing

Platform

CHAPTER 5

CRISPR sgRNA Genome Editing Platform




CRISPR-Cas9 is a major breakthrough after TALEN. The CRISPR-Cas9 system is a prokaryotic immune system driven by RNA that triggers resistance to foreign genetic elements and provides a form of acquired immunity. sgRNA serves as a guide to recognize and combine the target with DNA, allowing activated Cas9 nuclease to attach, cleave, and modify the PAM downstream target DNA. CRISPR-based genome editing technology has made it easier and faster than ever to alter specific DNA sequences in the genome or to perform genome-wide functional screening tests to identify genes involved in a particular phenotype. CRISPR-Cas9 is considered the most promising tool in gene modification. Synbio Technologies's CRISPR-Cas9 technology, combined with our patented synthesis and assembly platform, allows for efficient editing of genes/genomes.

Competitive Advantages


- ▶ Professional team offers customized design of sgRNA for gene family or genome region, without gene sequence or species restrictions
- ▶ Comprehensive services, including ready-to-use sgRNA, sgRNA vector construction, sgRNA library construction services, etc.
- ▶ One-stop solution from design, construction, to gene editing

CRISPR-Cas9 Applications



5.1 CRISPR sgRNA Design Center

The sgRNA design center of Synbio Technologies provides services such as single sequence sgRNA design or whole-genome sgRNA library design, downstream verification, and stable cell line construction. We offer a one-stop solution for CRISPR-Cas9 projects to achieve high genome editing efficiency.



Competitive Advantages

Professional Design Center

We provide single sgRNA, multiple sgRNA, and sgRNA library design services

Experienced

Well-designed sgRNA sequences with lowest off-target rate; Professional analysis is also provided to ensure the success of each project

Applicable in Many Species

Database covers sgRNA design for >20 different species; Specific sgRNA design is also available upon request for certain species

Custom Vector Construction



sgRNA can be designed in any vector that customers require

5.2 CRISPR System Activity Detection

In the CRISPR-Cas9 system, different guide RNA (gRNA) can be designed based on the target gene. The DNA cutting efficiency of Cas9 can vary depending on the gRNA and target gene. It is therefore important to be able to accurately verify the activity of the CRISPR-Cas9 system in a given assay, in order to maximize the cleavage and gene knock-in/knock-out efficiency of CRISPR-Cas9. Synbio Technologies provides SSA activity assays, *in vitro* cleavage activity assays, and endogenous activity assays to detect gRNA activity.


SSA Activity Assays

Single-strand annealing (SSA) assays can verify whether the target gRNA plasmid can mediate Cas9 to cut naked DNA. It is a common method to assess the activity of CRISPR-Cas9. The plasmid-based SSA assay contains two inactive luciferase-encoding fragments, which contain a stop codon and a segment of gRNA in the middle of the target sequence. If CRISPR-Cas9 can recognize and cleave the target site, the resultant double-stranded break can be repaired by the SSA mechanism, rescuing the now-active luciferase gene. Thus the SSA assay can be used to predict the cleavage activity of CRISPR-Cas9 by fluorescence detection.




Verification of *in vitro* Cleavage Assay

The target DNA sequence is digested with Cas9 *in vitro* to obtain two DNA fragments. Agarose gel electrophoresis allows analysis of gRNA-mediated DNA cleavage by comparing the yield of the smaller, digested fragments to the remaining undigested target DNA. Verification of *in vitro* cleavage is another way to determine the activity and efficiency of CRISPR-Cas9.



Endogenous Activity Assay

Mismatch endonuclease can also be used to detect sgRNA activity. The suitable primers are first designed on both sides of the target, followed by the PCR amplification of the DNA containing potentially mismatched mutation sites. T7E1 endonuclease can recognize and cleave mismatched DNA heteroduplexes. The digested product is finally analyzed by agarose gel electrophoresis to estimate the ratio of the cleavage and un-cleaved bands, reflecting the activity of CRISPR-Cas9.



Service Specifications

Services	Details	TAT	Deliverables
SSA activity assay	<ul style="list-style-type: none">gRNA target primer design and synthesispSSA-luc vector constructionCell transfection	Inquiry	Assay Report
<i>In vitro</i> cleavage activity assay	<ul style="list-style-type: none"><i>In vitro</i> transcription gRNACRISPR-Cas9 <i>in vitro</i> cleavage reaction		
Endogenous Activity Assay	<ul style="list-style-type: none">Genome extractionTarget primer designT7E1 enzyme digestion		

Comparison of different activity assay services

Assay Services	Differentiation
SSA activity assay	Easy to design, simple to use, cost-effective
<i>In vitro</i> cleavage activity assay	Easy to use (only target design is needed), cell-free assay system, cost-effective
Endogenous activity assay	Fast turnaround time, accelerated validation process

5.3 Ready-to-Use sgRNA Synthesis

CRISPR-Cas9 is one of the most convenient methods of gene editing technology, and has been widely used for editing genes for many different species. The traditional methods to deliver sgRNA into the cell are by plasmid transfection and lentiviral infection, which both have the risk of gene insertion and immune response. Currently, the novel way of directly transferring the ribonucleoprotein (RNP) formed by Cas9 and sgRNA into cells has the advantages of being safer and faster with a lower off-target effect and higher editing efficiency. This has gradually become a much more efficient way to use CRISPR/Cas9 technology.

To improve the efficiency of the CRISPR-Cas9 system, Synbio Technologies has developed ready-to-use sgRNA synthesis services. The ready-to-use sgRNA products can be directly transfected into cells or animals, avoiding the drawbacks of potentially non-degraded plasmids.

Chemical Synthesis of CRISPR-Cas9 sgRNA

Synbio Technologies uses a chemical synthesis strategy to provide our sgRNA synthesis service with sequences and desired modifications that are 100% correct. The synthesized sgRNA has more stability, low toxicity, and high editing efficiency. It's an excellent choice for your gene editing experiments.




Figure 1. Cas9 complex with sgRNA editing the genomic DNA




Figure 2. Secondary structure of the sgRNA's target and scaffold sequence

Competitive Advantages

- High Efficiency: The sgRNA is chemically modified, which will increase the stability and editing efficiency in vivo.
- Good Stability: Unrivaled control of sgRNA specifications ensures batch-to-batch consistency and traceability.
- Convenient: The sgRNA is ready-to-use and covers the sequence and function of both crRNA and tracrRNA, without annealing.

Service Specifications

Length (nt)	Order Size	Modification	Purification	Turnaround Time (Business Days)	Deliverables
95-105	1-10 OD	2'-OMe and Phosphorothioate	HPLC	7-10	1. Lyophilized RNA 2. Certificate of Analysis (COA) 3. QC Report (MS/HPLC)

In Vitro Transcription of CRISPR-Cas9 sgRNA

Synbio Technologies can also provide one-stop services including sgRNA target design, DNA template synthesis of the sgRNA, sgRNA *in vitro* transcription, and sgRNA purification, providing customers with sgRNA products that can be transferred into cells directly. According to publications, *in vitro* transcription of sgRNA has successfully edited the genes of many different species including zebrafish, mouse, filamentous fungi, etc.

Competitive Advantages

One-Stop Solution

Synbio Technologies provides integrated services from sgRNA target design to high purity ready-to-use sgRNA production

Fast Turnaround Time

In just 3 business days, Synbio Technologies will deliver up to 20 µg of customized ready-to-use sgRNA

Convenience

Ready-to-use sgRNA can be directly injected into animals or transfected into cells, improving the efficiency of gene editing experiments

Work Flow



Service Specifications

Services	Specifications	Deliverables
Ready-to-use sgRNA synthesis	sgRNA design DNA template synthesis <i>In vitro</i> sgRNA transcription and purification	sgRNA COA
Syno® negative control sgRNA	Negative control sgRNA <i>In vitro</i> sgRNA transcription and purification	Syno® negative control sgRNA1 Syno® negative control sgRNA2 Syno® negative control sgRNA3

Case Study

Synbio Technologies has designed several sgRNAs to target several genes in the mouse genome and have performed *in vitro* transfection with these constructs. The experimental period was shortened to 2 days and the sgRNA amount was increased to 10-20 µg. This change could potentially result in a significant jump in efficiency for synthetic biology experiments utilizing CRISPR-Cas9.



Result*
1. Clone DNA template into pUC57. The sequencing result (Fig.1) coincided with the designed sequence.




Fig.1 Comparison between blunt end ligation result and designed sequence

2. Agarose gel electrophoresis of sgRNA obtained by *in vitro* transcription. Clear bands shown in Fig.2.





Fig.2 Agarose gel electrophoresis of sgRNA

3. sgRNA Verification: Transcribe sgRNA into cDNA, design sgRNA amplification primer, and obtain the complementary DNA sequence by PCR reaction. Clone DNA sequence into pUC57 vector. The sequencing result (Fig.3) showed the sgRNA sequence is correct.



*The template of Lane 1 is reverse transcribed cDNA, the template of Lane 2 is sgRNA digested by DNase I, and the template of Lane 3 is *in vitro* transcribed DNA.

Fig.3 sgRNA sequence verification using agarose gel electrophoresis

5.4 CRISPR sgRNA Vector Construction

The CRISPR-Cas9 system is one of the most popular types of gene editing technologies and has been used in many research projects over the last few years. Due to its ease of use and strong applicability, CRISPR-Cas9 has been applied to biology fields including laboratory research, biopharmaceutical research, and genetic engineering of flora and fauna. This technology can also be applied to the construction of genetically modified animal models, construction of gene-regulated cell lines, breeding of plants, cancer research, and treatments of other serious diseases/screening of relevant drug targets. Synbio Technologies offers sgRNA plasmid construction services for several different animal and plant models. We guarantee fast and high-quality sgRNA plasmid construction, tailor-made to meet a large variety of requests and specifications set by our customers.

Competitive Advantages

Many Plasmid Options

We can provide several validated plasmids for different species to help you best fulfill the requirements of your experiment

High Success Rates

Several target sites can be designed for different genes which enhance the success rate of plasmid construction

Plasmid Validation Service


We also provide sgRNA activity validation services, which can enhance the success rate of experiments

Service Specifications

Processes	Details	
Confirm species and vector	Different species and different vectors	
Target site design	We usually design 3-5 target sites on common exon for different transcription products. The important domain should be damaged as far as possible which causes loss-of-function in the target gene.	
Synthesis	Should check the target site for single nucleotide polymorphisms (SNPs) before synthesis. (If there are SNPs, T7E1 enzyme cannot be used to detect mutations)	
Target cutting activity detection	SSA activity assays	<i>in vitro</i> cleavage activity assay
Endogenous active validation (Conditional choice)	If the target cell has a high transfection rate, like 293T cells, the genome can be extracted after 72 h. Afterwards, amplify target sequence, verify T7E1 enzyme digestion, and validate by sequencing	

5.5 ssDNA Synthesis

ssDNA (single-stranded DNA) is usually used as a homologous recombination repair (homology directed repair, HDR) template to achieve site-specific repair and long fragment knock-in in gene editing. Synbio Technologies provides high-purity, 100% accurate long ssDNA synthesis services.



Competitive Advantages

Sequence Guarantee

Sanger sequencing ensures 100% sequence accuracy

Non-toxic

The enzymatic approach ensures that ssDNA is non-toxic and has higher homology repair and knock-in efficiency

Highly Customized

Customised to produce the desired ssDNA for different application use


Service Specification

Length (nt)	Order Size (μg)	Turnaround Time (Business Days)	Deliverables
200-400	2-200	Start from 10-15	• Lyophilized ssDNA • Sequencing chromatogram • Certificate of analysis (COA)
400-1,000		Start from 10-15	
1,000-2,000		Start from 10-15	
3,000-4,000		Start from 15-20	
>4,000		Inquiry	

5.6 CRISPR Sequencing

When CRISPR technology is used for gene editing, a large number of target sites or samples are generated, which is suitable for high-throughput sequencing. That is, primers are designed near the target gene region, PCR amplification is performed, and then the amplicons are sequenced. Unlike CRISPR sequencing, traditional Sanger sequencing cannot quickly and efficiently identify the editing efficiency of a large number of genes. Synbio Technologies's amplification and sequencing technology can be used to detect CRISPR/Cas9 genome editing results and analyze off-target effects, helping researchers quickly locate gene mutation sites and study gene function.

CRISPR Amplicon Sequencing Process



Competitive Advantages

High Coverage

Each reaction can perform multiplex analysis of hundreds to thousands of amplicons

High Flexibility

It can be used for all kinds of mutation verification and screening genetic variation

Highly Customized

Analysis services can be customized

Cost-effective

Compared with whole-genome sequencing, it can reduce sequencing cost and turnaround time

Service Specification

Amplicon Size	Platform Configuration	Turnaround Time	Deliverables
70-280 bp	Illumina 2×150 bp	3-4 weeks	• Sequencing chromatogram • Data analysis report
280-480 bp	Illumina 2×250 bp	4-5 weeks	

CHAPTER 6

Protein Expression and Purification

With our experienced research team and leading protein expression platform, Synbio Technologies provides our customers a variety of protein expression and purification platforms, including bacterial, yeast, insect, and mammalian expression systems. By utilizing our powerful gene synthesis platform, Synbio Technologies can generate high-purity recombinant proteins ranging from milligrams to grams. Combined with our patented NG™ Codon Optimization Technology, the protein expression level can be significantly improved. Simply submit your gene or protein sequence and your high-quality purified recombinant protein will arrive to your bench in as soon as 3 weeks.

Competitive Advantages



► NG™ Codon Optimization Technology

Increase the protein expression level in most host cells

► Comprehensive Expression Systems

Bacteria, yeast, baculovirus/insect cell, and mammalian cell expression systems, plus various vectors for each system


► One-stop Services

Streamlined process from gene synthesis, codon optimization, to protein expression and purification


► Scalable Synthesis Capacity

Provide your desired proteins from milligram to gram scale

Service Procedure



Bacterial Expression System



Yeast Expression System



Baculoviral-Insect Cell Expression System



Mammalian Cell Expression System

Protein Expression and Purification

Protein Expression and Purification Services

6.1 NG™ Codon Optimization Technology

Codon optimization is a method of utilizing preferred codons to increase overall sequence stability and resulting protein expression. This is accomplished by avoiding rare codons with low utilization, simplifying the mRNA's secondary structure, and adjusting the overall GC content. The increased commercialization of gene synthesis has greatly promoted the application and influence of codon optimization. Synbio Technologies's NG™ Codon Optimization Technology optimizes complex sequences resulting in a significantly higher protein expression level.

Competitive Advantages

✓ Raise the protein expression level significantly

✓ Gene synthesis codon optimization is available in all hosts

Different organisms have their own preferred codons. For example, the prioritized termination codon of yeast and mammalian are UAA and UGA, respectively. Specific gene synthesis codon optimizations are applied for different organisms. Specific codon optimizations are provided in correlation to the hosts that clients provide.


✓ Excellent results even in the most difficult sequences

NG™ Codon optimization is capable of optimizing the most complex sequences such as hairpin sequences, highly repeated sequences, and extremely long sequences.

✓ Intellectual property (IP) protection

All of your IP-related rights are fully protected by Synbio Technologies. The DNA sequence and its related information are kept strictly confidential in the entire process. All materials and mid-products including synthesized oligos, genes, plasmids, and the cells will be destroyed in one month after delivery.

Case Study



T: Test
C: Control
M: Marker

Target proteins were reproducibly expressed in *E. Coli* system after codon optimization by NG™ Codon Optimization Technology.

6.2 Guaranteed Bacterial Protein Expression

With years of experience in protein expression, Synbio Technologies has developed a cost-effective bacterial expression platform to deliver high-quality proteins to our customers. By using our self-developed NG™ codon optimization technology and mature expression system, the expression level in bacteria can be significantly improved.

We offer guaranteed services, starting with our proprietary codon optimization, followed by gene synthesis, vector construction, all the way to protein expression and purification to ensure the expression of the target protein. If we fail to give you your desired protein, there will be no charge.

You only need to provide us with the sequence of the target protein. Then, we can deliver the specified quantity and purity of the target protein within a minimum of 4 weeks, ensuring your downstream experiments.

Competitive Advantages

01

Proprietary NG™ Codon Optimization Technology
Significantly improve protein expression level and solubility, with a success rate of up to 95% in the expression system of *E. coli*

02

Comprehensive Solutions
Provide one-stop solutions from gene synthesis to protein expression and purification

03

High Quality Guaranteed
There is no charge if we fail to deliver you desired proteins

04

Competitive Price and Fast Turnaround Time

Service Specifications

Services	Descriptions	Deliverables	TAT
Gene Synthesis	<ul style="list-style-type: none">Codon optimization and gene synthesisSubcloning into appropriate expression vectorVarious tags are available, such as His, SUMO,TRX,GST...		≤2 weeks
Protein Expression Testing and Purification, Tag Removal (Optional) & Refolding (Optional)	<ul style="list-style-type: none">Transform plasmids into appropriate bacterial expression strainProtein expression evaluationScale up protein expressionPurification to reach desired protein amounts and purityTag removal and separation of tag-free protein (if requested)Refolding in case of insoluble protein (performed only if required)	<ul style="list-style-type: none">2-5 µg lyophilized plasmid DNA1-5 mg protein productsPurity options: 85%, 90%, 95%Endotoxin level options: ≤0.1 EU/µgCertificate of analysis	≤3 weeks
QC & Delivery	<ul style="list-style-type: none">SDS-PAGEWestern Blot (if requested)Bradford assay for quantitationProtein delivery		≤1 weeks

* Suitable for proteins with molecular weight between 15-120 KD, excluding membrane proteins and toxic proteins.

* For gram level large-scale protein production, please contact quote@synbio-tech.com.

6.3 Yeast Protein Expression System

Prokaryotic expression is mature and easy to be operated, but when it comes to the proteins of eukaryotic organisms, error folding is a frequent occurrence since glycosylation and folding partners are involved. Fortunately, the expression system of *pichia pastoris* is complementary to this situation to a certain extent. The *Pichia pastoris* expression system is a simple and economic system among eukaryotic systems. It has many advantages in protein post-processing, folding, and post-translational modification. Synbio Technologies has abundant resources in vectors and host cells, and can effectively improve protein yield and provide you high-quality and considerate yeast expression protein services.

Service Specifications

Services	Descriptions	Deliverables	Turnaround Time
Gene Synthesis (Optional)	<ul style="list-style-type: none">Codon optimizationGene synthesisVector construction	<ul style="list-style-type: none">2-5 µg lyophilized plasmid DNACertificate of analysis (COA)	Starting from 10 BDs
Electrotransformation and Screening Positive Clones (PCR method)	Linearize the plasmid and transform into pichia pastoris cells to screen positive clones	Screening report	2-3 weeks
Pilot Protein Expression	Three positive clones will be used to verify the expression, then perform expression analysis validation (WB)	Expression evaluation report	4-6 weeks
Protein Expression and Purification	<ul style="list-style-type: none">1 L or larger scale expressionOne or more purifications stepsPurity will be detected by SDS-PAGE or Western blot (Optional)	<ul style="list-style-type: none">Purified protein productsCertificate of analysis (COA)	3-4 weeks

6.4 Baculovirus-Insect Cell Expression System

The baculovirus insect cell system offers several unique advantages over the *E. coli* expression system, such as post-translational modifications, higher yield for secreted proteins, and improved solubility. It is appropriate for expressing proteins that maybe harmful to mammalian cells, such as kinases and toxic proteins. With years of experience in protein expression, Synbio Technologies possesses abundant experience in baculovirus insect cell protein expression, particularly the expression of toxic proteins, kinases, large proteins, membrane proteins, and intracellular proteins.

Competitive Advantages

01

Proprietary NG™ Codon Optimization Technology

Significantly improves protein expression levels and success rate

02

Multiple Purification Strategies

Affinity chromatography, ion exchange, hydrophobic chromatography, gel filtration, etc., providing diversified choices of protein purities

03

Comprehensive Protein Identification Services

SDS-PAGE, Western Blot, HPLC, etc

Service Specifications

Services	Descriptions	Deliverables	Turnaround Time
Gene Synthesis (Optional)	<ul style="list-style-type: none">Codon optimizationGene synthesisVector construction	<ul style="list-style-type: none">2~5 µg lyophilized plasmid DNACertificate of analysis (COA)	Starting from 10 BDs
Baculovirus preparation	<ul style="list-style-type: none">Generation of recombinant Bacmid DNAGeneration of P1 stock, P2 stock, and determination of virus titer	Experiment report	3-5 weeks
Pilot Protein Expression and Purification	<ul style="list-style-type: none">P2 generation virus infects SF9 cellsExpression analysis validation (WB)Ni column affinity purification	<ul style="list-style-type: none">A small amount of proteinExpression evaluation report	2-4 weeks
Protein Expression and Purification	1 L or larger scale expression	<ul style="list-style-type: none">Purified protein productsCertificate of analysis (COA)	1-2 weeks

6.5 Mammalian Cell Expression System

Protein with post-translational modifications and appropriate protein folding are critical for many downstream applications. Synbio Technologies's mammalian expression platform offers a variety of services, including eukaryotic protein expression, recombinant antibody production, etc. Combined with our advanced protein purification facilities, we can deliver recombinant protein products with high purity $\geq 98\%$ and endotoxin level less than 0.001 EU/ug. With our Avitag technology, we also offer a variety of biotinylated proteins and antibodies to facilitate your assay development.

Competitive Advantages

High Efficiency

Various self-developed high expression vectors, meeting the different screening needs of customers

High Purity and Quality

Many of our antibody products have a purity higher than 95%, with comprehensive quality control

Low Endotoxin Level


Endotoxin level is less than 0.005 EU/ μ g

Service Specifications

Services	Descriptions	Deliverables	Turnaround Time
Gene Synthesis (Optional)	<ul style="list-style-type: none">Codon optimizationGene synthesisVector construction	<ul style="list-style-type: none">2-5 μg lyophilized plasmid DNACertificate of analysis (COA)	Starting from 10 BDs
Pilot Protein Expression and Purification	<ul style="list-style-type: none">Transient transfectionPilot protein expressionExpression analysis validation (WB)	Expression evaluation report	2-3 weeks
Protein Expression and Purification	<ul style="list-style-type: none">200 ml, 1 L or larger scale expressionPurity Identification (PAGE)	<ul style="list-style-type: none">Purified protein productsCertificate of analysis(COA)	2-3 weeks

Selection of Fluorescent Proteins, Promoters, and Selectable Markers


Different promoters have different activation activities in different cells. For example, in hematopoietic stem cells, SFFV has much better activation activity than CMV. There are also differences in tissue-specific expression. AAT is highly expressed in liver tissues and weakly expressed in other tissues. hSyn is highly expressed in nerve cells and weakly expressed in other tissues. Synbio Technologies provides a variety of fluorescent proteins, promoters, and mammalian expression vectors with selectable markers to meet the needs of cell-specific screening.



6.6 Peptide Synthesis

With our experienced research team and leading synthesis platform, Synbio Technologies provides our customers with a full spectrum of peptide services ranging from custom peptide synthesis to peptide library synthesis. Our custom peptides and peptide libraries are available with flexible modification, conjugation, and isotopic labeling options. Up to date, our success rate of peptide synthesis has achieved more than 98%. The high quality peptide synthesis products provided by Synbio Technologies can be applied to drug development, protein functional analysis, monoclonal antibody preparation, antibody-antigen interaction, enzyme specificity research, etc.

Competitive Advantages



Service Specifications

Services	Turnaround Time	Deliverables
Custom Peptide Synthesis	2-3 weeks, 1-2 weeks for urgent orders	<ul style="list-style-type: none">• Peptide Products• HPLC Reports• MS Reports• Certificate of Analysis
Peptide Modification		
Peptide Library		

* Please contact quote@synbio-tech.com for more information.

Types of Peptide Libraries

Overlapping Peptide Library

DRVYIHPFHLWIEG
DRVYIHPFHL
RVYIHPFHLWIE
VYIHPFHLWIE
VYIHPFHLW
VYIHPFHL

Truncation Peptide Library

DRVYIHPFHLWIEG
DRVYIHPFHLW
DRVYIHPFHL
DRVYIHPFH
DRVYIHPF
DRVYIHP

Alanine Peptide Scanning Library

DRVYIHPFHLWIEG
DRVYIHPFHLWIEG
DRVYIHPFHLWIEG
DRVYIHPFHLWIEG
DRVYIHPFHLWIEG
DRVYIHPFHLWIEG

Random Peptide Library

DRVYIHPFHLWIEG
DRVYIHPXXLWIEG
DRVYIHPXXHLWIEG
DRVYIHPFHXXEG
DRVYIHPFXXWIEG
DRVYIHPXXLWIEG

Scrambled Peptide Library

DRVYIHPFHLWIEG
DVPYIHRFGEWILH
DRGWHPFHLWIEG
HLVPIHYFDRWIEG
RIVYDHPFHLWIEG
DHFYIHPFVRWIEG

Positional Peptide Library

DRVYIHPFHLWIEG
DRVYIHPFHLWIEG
DRVYIHPFHLWIEG
DRVYIHPFHLWIEG
DRVYIHPFHLWIEG
DRVYIHPFHLWIEG

Individual peptides can be divided into several fragments that overlap. The resulting overlapping peptide libraries can then be used for processes including continuous and linear epitope mapping.

The library of truncated peptides can be used to predict the minimum length amino acid that could be used to achieve optimal epitope activity.

Each amino acid is substituted individually and systematically for alanine. This library enables quick determination of each individual amino acid's contribution to the peptide's functionality.

Selected positions are substituted with all 20 natural amino acids simultaneously, which might increase peptide activity.

Scrambled libraries have the highest variation of any peptide library. The resulting peptides are used generally as negative controls to show that a specific sequence is critical to the protein function or activity. It is also a random screening tool used to find new leads.

Selected sites or regions within a peptide sequence are replaced systematically with all other natural amino acids. This allows for amino acids that enhance peptide activity to be identified.

Antibody Services

CHAPTER 7

Antibody Services

Antibody engineering refers to the modification and re-assembly of antibody genes through recombinant DNA and protein engineering technologies. The engineered genes are then further transfected into receptor cells and expressed accordingly. Additionally, antibody engineering also includes the reformatting of known antibodies via cell fusion and chemical modification. Both antibody engineering approaches yield new artificial antibody molecules that retain or newly add the specificity and biological activity of natural antibodies, or eliminate unrelated structures. For this reason, antibody engineering has great potential to further the capabilities of pharmaceutical research and development.

Synbio Technologies provides a broad range of DNA technologies to antibody services. Our technical platforms span all stages of antibody discovery by integrating “DNA Reading” (sequencing & discovery), “DNA Writing” (design & synthesis), “DNA Editing” (optimization & humanization). In addition, we can provide high-quality and cost-effective recombinant antibody products to support the development of biobased medicines in academia as well as in the biotechnology and pharmaceutical industries.

Competitive Advantages


Syno® Ab antibody design platform uses proprietary molecular simulation algorithms to reduce the time and cost of antibody discovery and shorten the development cycle of antibody drugs

Advanced Syno® synthesis technology platforms provide customers with fast and efficient recombinant antibody production

Geno™ Ab antibody discovery platform and our professional bioinformatics teams provide our customers with highly accurate antibody sequencing and fast antibody discovery

Customized service to meet your personalized needs


Provide one-stop solutions from sequencing and synthesis to antibody preparation



7.1 Monoclonal/Hybridoma Antibody Sequencing

Obtaining the hybridoma monoclonal antibody (mAb) sequence is one of key requirements of recombinant antibody engineering. Monoclonal antibody sequencing includes hybridoma cell line sequencing and clonal B cell sequencing. Both of these kinds of cells can produce a single antibody against a unique antigen after screening. Synbio Technologies uses 5'RACE and Sanger sequencing technologies to provide both accurate and efficient hybridoma/clonal B cell sequencing services across a variety of species.

Competitive Advantages



Hybridoma Sequencing Process

mRNA isolated from hybridoma samples is subjected to RT-PCR with unique species-specific primer sets to amplify the target regions. Synbio Technologies's antibody variable region sequencing service has identified more than 10,000 hybridoma cell lines. More than 95% of the samples have been successfully sequenced with both the heavy chain and light chain identified. Even with technically difficult samples, Synbio Technologies's approach has demonstrated a high success rate in numerous cases. Synbio Technologies's services, coupled with our online data tracking and management system, provides customers with absolute confidence at every step.

CHAPTER 7

Antibody Services


Standard Hybridoma Sequencing Process



Service Specifications

Services	Descriptions	Turnaround Time	Deliverables
Standard mAb Sequencing	Standardized work flow for hybridoma sequencing	2-3 weeks	<ul style="list-style-type: none">• Raw sequence data (100% accuracy guarantee)• Final sequence report with variable heavy chain and light chain sequences or full-length sequence and specific regions annotated• Plasmid contains antibody coding sequence (upon your request)
Express mAb Sequencing	Rapidly identify antibody sequences from the target cell lines	1 week	

* For antibody sequencing services: customers will need to submit hybridoma cell samples (>1×10⁶ cells) and indicate the subtype of your monoclonal antibody.



7.2 Immune Repertoire Sequencing

"Immune Repertoire" refers to the set of all functionally diversified B cells and T cells in the circulatory system of an individual at any given time. Immune repertoire sequencing (Immuno-Seq) targets T/B lymphocytes and utilizes 5' Rapid Amplification of cDNA Ends (RACE) Technology. 5' RACE technology amplifies complementarity-determining regions (CDRs) of B cell receptors (BCRs) or T cell receptors (TCRs). The CDRs are targeted because they control the diversity of the corresponding BCRs and TCRs. Combined with high-throughput sequencing technology, Immuno-Seq provides a comprehensive evaluation of immune system diversity.

Synbio Technologies provides our customers with a complete antibody repertoire sequencing system including immune repertoire sequencing, antibody CDR region sequencing, etc. Our Immuno-Seq System is designed to assist researchers in the observation and analysis of both T and B cells with unprecedented specificity. The specificity allows the Immuno-Seq technology to be applied to disease surveillance, antibody production, vaccine research, health examination, and other related areas.

Competitive Advantages



Standard IR-Seq Process




CHAPTER 7

Antibody Services

Service Specifications

Service	Turnaround Time	Deliverables
Immune Repertoire Sequencing	6-8 weeks	<ul style="list-style-type: none">• NGS sequencing data• Bioinformatics analysis report

* For Immuno-Seq service: customers will need to submit one of the following: 1-2ml peripheral blood, splenocyte, selected immune cells, or extracted high quality RNA.



CASE STUDY

Synbio Technologies’s immune repertoire sequencing and bioinformatics analysis:

1. Library Construction

Synbio Technologies performed 5’ RACE reverse transcription of the extracted RNA, then amplified and constructed the library.

2. Library Quality Control


The Qubit 2.0 Fluorometer was used to perform preliminary quantification after the library construction. After dilution, insert size detection was performed on the library based on the Agilent 2100 bioanalyzer. Only when meeting expected requirements, the effective concentration of the library (>2nM) was accurately quantified by qPCR to ensure the quality of the library.

3. Library Sequencing

Different libraries were pooled and sequenced by Illumina Miseq/HiSeq according to the effective concentration and the required amount of target data. Four fluorescent labeled dNTPs, DNA polymerase, and adaptors were added to the flow cell for the amplification. When each sequence cluster extended and generated the complementary chain, the corresponding fluorescent signals were released and captured after each new dNTP was added. The optical signals were then transformed into sequencing peaks by professional software to obtain the sequence information of the fragments.

4. Data Analysis

Quality assessment of sequencing data.




7.3 De Novo Antibody Design and Production


Biological drugs usually include genetic engineering drugs, antibody engineering drugs, blood product drugs, vaccines, etc. Among these, antibody drugs are one of the fastest growing therapeutic drugs for humans. This increase in research and development has proved to be effective for treatment of many important diseases such as cancer, chronic inflammatory diseases, cardiovascular diseases, immune system diseases, and infectious diseases. Antibody research and development also plays an important role in medicine and diagnostics.

Synbio Technologies's proprietary antibody platform utilizes molecular simulation technologies to provide *de novo* antibody design and preparation services. We can mimic the biological behavior of antigenic molecules at the molecular level. We analyze antibody structures as a starting point to solve a series of problems, such as the theoretical design of an antibody. The combination of molecular simulation technology and our experimental approach can help our customers effectively reduce the total cost of antibody R&D while shortening their R&D cycle.

Competitive Advantages



De Novo Antibody Discovery by Syno® Ab



Service Specifications

Service	Details	Turnaround Time	Deliverables
Syno® Ab Antibody Design and Production	<ul style="list-style-type: none">• CDR construction• Antibody modeling• Epitope docking• Screening optimization	2-3 weeks	<ul style="list-style-type: none">• Antibody sequence report• Structure analysis of specific targets (antigens)• Molecular simulation and molecular dynamics ordering• Computational ordering of antibody library (human or specific species)• Leads with the affinity greater than 10⁻⁹ M• Antibody library or the synthesized leads antibodies

* Customers can provide the antigen name or gene ID directly. If you have the antigen's crystal structure, please contact us directly.


7.4 Antibody Humanization Service

Immunogenicity or a Human Anti-Mouse Antibody (HAMA) response is what happens when antibodies derived from xenogeneic sources are administered to humans. This can potentially decrease the effectiveness of immune-treatment through rapid clearance of the therapeutic antibody or immune-conjugate. In some cases, the patient may even develop an allergic sensitivity and be at risk of anaphylactic shock upon any future treatment containing these “foreign” antibodies. Due to this potential risk of immunogenicity and HAMA response, the humanization of antibodies from xenogeneic sources to become safer and more reliable humanized antibody drugs is critical.

Synbio Technologies has developed an integrated antibody humanization and affinity maturation platform based on advanced phage/yeast display technologies and progressive computer modeling and simulation techniques. Collectively, we have combined methods from bioinformatics and structural biology to optimize our approach toward integrated antibody humanization. Our antibody humanization services include *in silico* CDR-grafting, computer-aided design and optimization, high-throughput screening of phage-display antibodies, and Biacore surface resonance analysis. With this approach, we are able to provide humanized antibodies with diminished immunogenicity, which have comparable or better specificity, bioactivity, thermostability, and productivity than parental antibodies. Our affinity maturation service assures optimal humanized antibody output with an increased binding affinity by at least five-fold.

Featured Advantages

- ✓ **Extensive Experience:** Our team has practical experience with humanized antibody design, screening, and production to ensure the success of our customers' projects
- ✓ **In Silico Design:** Effective and reliable antibody modeling and simulating from multiple antibody databases
- ✓ **High-Throughput Screening:** Phage/yeast display platform designed to enhance the valid antibody yield
- ✓ **Guaranteed Affinity:** After affinity maturation, the final construct's affinity is improved by at least five-fold




Service Specification

Details	Deliverables	Turnaround Time
<ul style="list-style-type: none">Chimeric antibody constructionAntigen-antibody binding analysis	<ul style="list-style-type: none">Top 5 humanized antibody variants and corresponding DNA sequences.At least one variant with binding affinity comparable to parental antibody.Full report (PDF).	4-5 months
<ul style="list-style-type: none">Antibody modeling library constructionPhage/Yeast display antibody library construction		
<ul style="list-style-type: none">High-throughput antibody screening		
<ul style="list-style-type: none">Humanized antibody characterization and production		


7.5 Antibody Affinity Maturation Service

The binding affinity of an antibody is one of the key parameters to determine the efficacy and dosage of an antibody drug. Currently, the methods for improving the binding affinity of antibodies mainly focus on site-mutation, CDR rearrangement, DNA recombination, etc. With our unique and effective affinity maturation technology, we can quickly spot the key amino acid sites in antibody domains. This, in turn, increases the binding affinity of the antibody, reaching at least nM level.

Competitive Advantages



Phage-Display Screening Process




Service Specification

Details	Deliverables	Turnaround Time
DNA Shuffling	1. Antibody sequence and sequencing report 2. Top 1-3 of the optimized antibodies 3. Full analysis report (PDF)	6-7 months
Library Construction and Screening		
Antibody Production and Characterization		


7.6 Recombinant Antibody Production

With the development of antibody engineering technology, various engineered antibodies have been designed for cancer diagnosis and immunotherapy. Based on years of research experience and leading professional technologies in the field of antibody engineering, Synbio Technologies can provide a series of recombinant antibody preparation services, including single chain antibodies (scFv), antigen binding fragments (Fab), bispecific antibodies, etc.

Competitive Advantages



Recombinant Antibody Production Process




Service Specifications

Services	Service Details	Deliverables
Antibody Discovery	Hybridoma sequencing or immune repertoire sequencing (optional)	Antibody gene sequencing report
Gene Synthesis	<ul style="list-style-type: none">Gene synthesis and vector construction: IgG, Fab, scFv, VHH, BsAb, etcPlasmid preparation	2-5 µg lyophilized plasmid
Recombinant Antibody Expression and Purification	<ul style="list-style-type: none">Cell transfection & culture30 ml, 100 ml, 200 ml or larger scale recombinant antibody expressionAntibody affinity purification and purity detection (SDS-PAGE, SEC-HPLC), purity > 95%Endotoxin: <1 EU/mg	<ul style="list-style-type: none">Antibody productsQC report

7.7 Monoclonal Antibody Preparation

Hybridoma technology features the effective usage of innate functions of both immune and cancer cells to generate monoclonal antibodies (mAbs). Resulting hybridoma cell lines are produced to combat specific antigens of interest. Such technology has forged remarkable new approaches toward therapeutic mAb discovery, as well as disease diagnosis and prevention. Mice are the main host animals for monoclonal antibody production. Synbio technologies can provide customers with a full range of customized monoclonal antibody production services, including immunization, hybridoma development, antibody preparation, antibody function screening analysis, etc.

Competitive Advantages



Service Specifications

Milestones	Details	Deliverables	Turnaround Time
Animal Immunization	<ul style="list-style-type: none">Blood collection before immunizationImmunize 4 timesFinal blood collection	<ul style="list-style-type: none">Experimental report (including antiserum titer test results)The residual antigen	6-8 weeks
Cell Fusion & Screening	<ul style="list-style-type: none">Cell fusionSubcloning and cell amplificationCell screeningClone selection and freezing	<ul style="list-style-type: none">Experimental reports (including screening data)	5-6 weeks
Antibody Purification	<ul style="list-style-type: none">One hybridoma cell cultureAntigen affinity purification or Protein A/G affinity purificationAntibody purity test (SDS-PAGE)Antibody titer detection (ELISA, titer > 1:100000)	<ul style="list-style-type: none">1-3 mg monoclonal antibody products, titer >1*10⁵Full experimental report	4-5 weeks

7.8 Polyclonal Antibody Preparation

Polyclonal antibodies are a mixture of antibodies that recognize different epitopes of a specific antigen and can be produced by collecting antiserums after artificially immunizing the host. Because it can recognize multiple epitopes, it has a higher affinity and sensitivity than monoclonal antibodies. Polyclonal antibodies are also more tolerant to small antigen changes. With rich experience in antigen design and mature antigen/antibody preparation technology, Synbio Technologies can provide one-stop customized services from antigen design, gene synthesis, antigen preparation, animal immunity, to high throughput antibody production. In addition, some host animals (such as mice, rabbits, goats, etc.) and various kinds of antigens (such as peptide antigens, protein antigens, small molecule antigens, etc.) are available. Our technical team will evaluate your requirements and tailor your project to meet your needs. Ultimately, we can deliver high-quality polyclonal antibodies for your scientific research.

Competitive Advantages



Service Specifications

Milestones	Details	Deliverables	Turnaround Time
Antigen Preparation	Antigen provided by the customer: 1. The peptide 10 mg 2. Recombinant protein 5-8 mg, more than 85% purity 3. Small molecular antigens	The experimental report of antigen preparation	2-4 weeks
	Synthesized or prepared by Synbio Technologies: 1. Design and synthesize polypeptide, coupling with KLH/BSA/OVA 2. Codon optimization, gene synthesis, recombinant protein expression and purification		
Animal Immunization	Immunize 3-4 times	/	7-8 weeks
Antibody Purification	Antigen affinity purification or Protein A/G affinity purification	0.5-10 mg antibody products, purity>80%	1-2 weeks
Quality Control	<ul style="list-style-type: none">Antibody purity test (SDS-PAGE)Specific detection of antigen and antibody (Western Blot analysis)Antibody titer detection (ELISA, titer > 1:100000)	Experimental report	1-2 weeks

7.9 Single Domain Antibody (sdAb) Service

Complex structures and high production costs associated with monoclonal antibodies (mAbs) have limited their achievable clinical applications. For this reason, scientists focusing on antibody drug discovery have shifted their interests toward recombinant antibody fragments.

Recombinant antibodies, such as single domain antibodies (sdAb), are a new generation of engineered antibody fragment that consists of a single monomeric variable antibody domain. sdAbs are derived from camelid heavy chain antibodies which lack the light chain and CH domain of the heavy chain in the conventional Fab region. Therefore, sdAbs are smaller, more soluble, and have higher stability than conventional antibodies. sdAbs also feature higher permeability, lower immunogenicity, while being both simple and cost-effective. These advantages have made it possible for sdAbs to have a broader range of applications within the biotechnology, pharmaceutical, and immuno-therapeutic industries. Based on Synbio Technologies's extensive expertise in sdAb and advanced sdAb discovery technology, we can provide our customers with a comprehensive one-stop solution from antigen production to specific recombinant sdAb expression.

Competitive Advantages

01

Experienced antibody production team

02

Proprietary Universal Antigen™ immunization technology

03

Robust phage-display technology platform

04

Advanced yeast-display technology platform
(Large library capacity: ~ 10⁹)



05

Cost-effective and milestone payment plan

Service Specification

Milestones	Deliverables	Turnaround Time
Camel immunization (eg. dromedaries, camels, llamas, and alpacas)	1. 100-500 µg sdAb products 2. Experimental report	10-11 months
PBMC isolation & immune library construction		
Screening of phage-displayed library		
sdAb expression by eukaryotic or prokaryotic expression system		
Functional assay for sdAb		

Standard Process



CHAPTER 8

RNA Synthesis

RNA has been widely used in gene function analysis and novel therapeutic research/development, along with various other biological research fields. As the interest in RNA research has become more popular, the needs for RNA synthesis increase greatly.

Synbio Technologies has decades of experience in RNA research and development, our scientists can use both a chemical synthesis approach and IVT synthesis strategy of RNA to meet the unique requirements of every client. Synbio Technologies's RNA synthesis services include RNA oligo synthesis with various modifications and delivery specifications available, functional siRNA synthesis, miRNA synthesis, and *in vitro* transcription RNA synthesis. All our synthesized RNA sequences go through a standardized manufacturing process and undergo strict QC review. We promise to provide the most economical price to meet scientists' diversified needs without sacrificing on overall quality.

Competitive Advantages

01

Excellent quality and stability

02

Flexible synthesis scales and various modification types to meet your customized requirements

03

Competitive prices and short turnaround time

Applications

Synthesized RNA can be used in a number of applications, such as:

RNA probe synthesis for hybridization

dsRNA/siRNA synthesis for RNA interference


miRNA synthesis for gene function research

gRNA synthesis for gene editing

Synthesis of mRNA / cRNA for electrophysiological research

Synthesis of various types of RNA for noncoding RNA research

RNA virus genome synthesis for virus research



RNA Synthesis

8.1 Custom RNA Oligos

When the fragment of mRNA you need is small, between 5 bases to 120 bases, we recommend using a chemical synthesis strategy, which is a quick and more economical method to produce a high yield of chemically synthesized RNA for various applications.

Synbio Technologies provides high-quality and diverse chemosynthetic RNAs, including customizable single/double stranded RNA oligos, modified RNA oligos, stable siRNA, miRNA, long RNA, etc. Chemically synthesized RNAs can be synthesized to contain nearly any sequence and chemical modification desired, with up to gram-quantity yields possible. All of our RNA products' quality will be confirmed by MS and HPLC purity tests to ensure the high-quality output.

Competitive Advantages

01

Various types of customizable RNA synthesis services capable of providing RNA fragments up to gram-scale

02

Strict purification procedures of each batch of synthetic products by HPLC in an RNase-free environment, avoiding RNase contamination throughout the process, and a free HPLC detection report is provided

03

Professional technical team provides you with comprehensive RNA solutions

04

Commitment to keeping the details of the project strictly confidential, protecting your intellectual property rights

Service Specifications

Single-stranded RNA				
Length	Order Size	Purification	Turnaround Time (Business Days)	Deliverables
5-65 bases (General Sequences)	2-10 OD	HPLC	5-7	1. Lyophilized RNA 2. Certificate of Analysis (COA) 3. QC Report (MS/HPLC)
	25-100 OD	HPLC	6-8	
	≥250 OD	HPLC	Inquiry	
65-120 bases	Inquiry	HPLC	Inquiry	
Double-stranded RNA				
Each single strand is purified by HPLC and then annealed				

8.2 RNA Modifications

The modification of an RNA sequence can help the RNA to be delivered to the cell while increasing the stability of the RNA sequence, avoiding degradation by nuclease, such as phosphorothioate, 2'-OMe, and 2'-F modifications. The advantage of fluorescent labeling in RNA is to let the synthesized RNA be observed or located by fluorescence microscopy and confocal laser microscopy, while also using it as a nucleic acid detection tool to perform virus detection, disease diagnosis, and early tumor screening.

Synbio Technologies provides extensive modification options for synthesized RNA, making it convenient for customers to choose flexible and diverse modifications and fluorescent labeling according to the function of synthesized RNA sequences.

RNA Modifications and Fluorescent Labels

RNA Modifications		RNA Fluorescence Labels/ Quenchers Groups	Dual/Multiple Labels
Modified Bases	<ul style="list-style-type: none">• 2'-O-Methyl• 2' Fluoro• m6A	<ul style="list-style-type: none">• 5' Cy3• 3' Cy3• 5' Cy5• 3' Cy5• 5' FAM• 3' FAM• 5' ROX• 3' ROX• 5'HEX• 5'TAMRA• 5'TET• 3'BHQ-1• 3'DABCYL• 3'TAMRA	<ul style="list-style-type: none">• 5' FAM-3' BHQ1• 5' HEX-3' BHQ1• 5' TET-3' BHQ-1• 5' FAM-3' MGB• 5' FAM-3' TAMRA• 5' HEX-3' TAMRA• 5' TET-3' TAMRA• 5' HEX-3' MGB• 5' TET-3' MGB• 5' FAM-3' DABCYL• 5' HEX-3' DABCYL• 5' Phos-3' FAM• 5'TAMRA-3'Phos• 5'C6-Biotin-3'Cy5
Phosphorothioates	<ul style="list-style-type: none">• P-S		
Phosphorylation	<ul style="list-style-type: none">• 5'- PHO• 3'- PHO		
Amino Linkers	<ul style="list-style-type: none">• 5'-NH2 C6• 3'-NH2 C7		
Thiol	<ul style="list-style-type: none">• 5'-SH C6• 3'-SH C3		
Biotin	<ul style="list-style-type: none">• 5' Biotin		
Digoxigenin	<ul style="list-style-type: none">• 5' DIG• 3' DIG		

8.3 siRNA Synthesis

RNA interference (RNAi) using small interfering RNA (siRNA) is the simplest and most efficient way to knock down gene expression to study protein function in a wide range of cell types. It is also the best method for clinical treatment and drug development of RNAi.

Synbio Technologies provides double-stranded ordinary siRNA as well as modified or fluorescence-labeled siRNA. Available chemical modifications include phosphate backbone modifications, ribose modifications, and base modifications, which can increase the stability and duration of siRNA in vivo, or enhance base-base interactions to facilitate mRNA targeting.

Service Specifications

Products	Order Size (OD)	Purifications	Turnaround Time (Business Days)
NC FAM siRNA	1-2,500	HPLC/OPC	Starting from 3 business days
NC siRNA			
siRNA			
Modified or fluorescence labeled siRNA (2-OMe, Thiol, Chol, Biotin, FAM, Cy5, Cy3, etc.)	Inquiry		Inquiry

*For more information about chemical modification and fluorescent labeling, please contact us at quote@synbio-tech.com.

Guaranteed siRNA Service

Synbio Technologies also provides a guaranteed siRNA service. The three designed and delivered siRNA ensure at least one siRNA can effectively inhibit the expression of the target gene, and the inhibition efficiency will be more than 70% (The transfection efficiency will be at least 90%).

Products	Order Size (OD)	Purifications	Turnaround Time (Business Days)
Three siRNAs	2	HPLC/OPC	6-9
NC FAM siRNA	1		
NC siRNA	1		
Positive siRNA control	1		

8.4 miRNA Synthesis

miRNAs are single-stranded, small, non-coding RNA of 22 nucleotides that play important roles as endogenous gene regulators by mediating translation repression or promoting degradation of target mRNA. While the normal expression and function of miRNAs are vital for physiological processes, aberrant expression of miRNAs has been proved to be closely related to the occurrence of various cancers. Perturbation of endogeneous miRNAs help to study the function of specific miRNAs and hold the potential for therapeutics.

Synbio Technologies offers a wide range of high-quality miRNA synthesis products, including miRNA mimics/inhibitors, miRNA agomirs/antagomirs and miRNA negative controls, etc., to support your miRNA functional research.

Service Specifications

Products	Order Size (OD)	Purifications	Turnaround Time (Business Days)
NC FAM miRNA mimics	1-2,500	HPLC/OPC	Starting from 3 business days
NC miRNA mimics			
miRNA mimics			
NC FAM miRNA inhibitors			
NC miRNA inhibitors			
miRNA inhibitors			
miRNA agomirs			
miRNA antagomirs			

8.5 RNA *In Vitro* Transcription Service

Another strategy of RNA synthesis is by *in vitro* transcription, especially suitable for long RNA synthesis. *In vitro* transcription takes linear DNA sequences as a template, using T7, T3, or SP6 RNA polymerase to synthesize RNA from the DNA sequence. Synbio Technologies provides highly efficient and cost-effective *in vitro* transcription RNA synthesis services. You can provide us with the DNA templates, such as the plasmid and PCR product, or we can help to do the *de novo* gene synthesis of the target DNA sequences. In our design, we synthesize the T7 promoter before the target sequence, which is used to initiate the transcription reaction during the RNA synthesis process. For mRNA transcription, Synbio Technologies provides optional mRNA modifications, including cap structures, poly(A) tails, and various modified nucleosides.

Competitive Advantages

Sequence Optimization

Improve mRNA translation efficiency through codon optimization

Fully Customized

5 'cap, 3' polyA tail, and various nucleoside modification are available to enhance the stability and translation efficiency of mRNA and reduce its immunogenicity


High Performance

Standard experimental procedures, professional staff, and abundant expertise offer the most optimized and tailored solutions

One-Stop Service

Provide one-stop services from DNA template synthesis to RNA *in vitro* transcription

RNA *In Vitro* Transcription Process



Service Specifications

Service Types	Details		Purification Methods	Length	Turnaround Time	Deliverables
General RNA Synthesis	No UTR region, cap structure, poly(A) tail and nucleoside modification		<ul style="list-style-type: none">• Column Purification• LiCl Purification• Phenol/Chloroform Extraction	<9 kb	1-2 weeks	<ul style="list-style-type: none">• Liquid RNA products• COA files
mRNA Synthesis	Cap1 analogs + transcription 120 poly(A) tail using template	N1-Me-Pseudo UTP Ψ-UTP 5m-CTP AF488/Cy3/Cy5			1-2 weeks	
	Cap1 analogs + post-transcription with poly(A) tail					
	Enzymatic capped + transcription 120 poly(A) tail using template					
	Enzymatic capped + post-transcription with poly(A) tail					

CHAPTER 9

COVID-19 DNA Solutions

Since the outbreak of COVID-19, Synbio Technologies has set up a special coronavirus emergency project team, working against time to produce virus-related probes, genes, and protein products. We are currently providing comprehensive COVID-19 DNA solutions and working with our clients and partners from hospitals, scientific research institutes, and enterprises to fight the disease.

► Diagnostic Strategy

Probe/Primer Synthesis for Virus Detection

Synbio Technologies has prepared a series of primers and probes for COVID-19 detection in batches according to WHO and the CDC's latest official documents. We provide a free trial package and guarantee the supply of a million copies of nucleic acid diagnostic probes daily to speed up disease detection and research.

► Therapeutic Strategy

From Gene Synthesis to Drug Discovery

Our exclusive CI system can facilitate the accurate, fast, and high-quality gene synthesis of any sequence. Syno® Ab is our *de novo* antibody design & production platform based on the combination of DNA manufacturing with computer-simulated antibody molecule design. Syno® Ab helps researchers effectively reduce costs and time on designated target design, neutralizing antibody generation, and drug discovery.

► Preventive Strategy

Vaccine Development and Production

Syno® GPS is our advanced biotechnology transformation and application platform. We effectively serve the development and production of genetically engineered vaccines, nucleic acid vaccines, and therapeutic vaccines. Currently, we have completed the synthesis of all structural and non-structural protein genes / fragments of COVID-19. This facilitates the acquisition of viral recombinant proteins and vaccine development research.

Synbio Technologies has been using our efforts to realize the company's vision of "Genes for Life". We hope that, through our one-stop DNA solutions to COVID-19 and powerful synthetic technology platforms, we can accelerate virus research and help benefit each patient's life as soon as possible.



9.1 Probe Synthesis for COVID-19 Detection

In early January 2020, Synbio Technologies designed and synthesized detection probes for COVID-19 based on our independently developed Syno® qPCR primer/probe design and synthesis platform. In less than two weeks, we completed the synthesis of detection reagent materials, such as amplification primers, probes, and quality controls.

Competitive Advantages

01

Speed Guarantee
All COVID-19 related orders will be given priority to green channel production, with 24/7 professional online support

02

Quality Assurance
ISO 13485 standards of manufacturing, high-quality raw materials, zero cross-contamination, and perfect track records

03

Production Capacity Guarantee
Experienced synthesis experts, first-in-class instruments, more than one million probe kit materials supplied per day

Probe/Primer Synthesis for Virus Detection

According to the official website documents of WHO and the CDC, Synbio Technologies has prepared a series of primers and probe materials for COVID-19 detection in batches. We can provide a free trial package and guarantee the supply of a million copies of nucleic acid diagnostic probes daily to speed up disease detection and research.




Figure 1. Relative positions of amplicon targets on SARS-CoV an 2019-nCoV genome.
ORF: Open reading frame
RdRp: RNA-dependent RNA polymerase. Numbers below amplicon are genome positions according to SARS-CoV, NC_004718.

► SARS-CoV-2 Real Time RT-PCR Detection Assay

To help facilitate COVID-19 research and contribute to fight the virus, Synbio Technologies has manufactured a detection assay based on the CDC's Real-Time RT-PCR protocol for detection of the 2019 novel coronavirus. This assay has successfully detected SARS-CoV-2 in positive control samples N gene and human RNase P (RNP) gene to assess the specimen quality. It should be used for research use only (RUO).

Probe Assay Specifications

Box #1: Primers and Probes

Reagent Label	Part #	Description	Quantity /Tube	Reactions /Tube
SARS-CoV-2_N1 SARS-CoV-2_N2 SARS-CoV-2_N3	RV202001 RV202002 RV202003	SARS-CoV-2_N1/N2/N3 Combined Primer/Probe Mix	0.56 nmol	24
RP	RV202004	Human RNase P Forward Primer/Probe Mix	0.56 nmol	24

Box #2: Positive Control

Reagent Label	Part #	Description	Quantity	Note
RNAPC	RV202005	SARS-CoV-2 Positive Control (RNAPC)	1 tube	Providing (24) 5 µL test reactions

Features

- ✓ Based on CDC protocol, ISO 13485 standards of manufacturing reagents.
- ✓ TaqMan™ Real-Time RT-PCR assay primers and FAM probes for higher specificity.
- ✓ Targeting three targeted sequences of N gene in GenBank sequences NC_045512.2.
- ✓ Positive control included already, no extra preparation required.

STATEMENTS

- This assay is for Research Use Only (RUO) and has not been tested on clinical samples.
- This assay is designed based on CDC protocol and public SARS-CoV-2 genome sequences.
- Different sample extraction methods, sample quality, operation parameters, data processing, and other un-tested factors might impact this assay's performance.

CHAPTER 9

COVID-19 DNA Solutions



► Primer/Probe Set 01

Name	Description		Oligonucleotide Sequence (5'>3')	Length	5'mod	3'mod
RdRP gene	1	Primer RdRP_SARSr-F2	GTGARATGGTCATGTGTGGCGG	22		
	2	Primer RdRP_SARSr-R1	CARATGTAAASACACTATTAGCATA	26		
	3	Probe RdRP_SARSr-P2	CAGGTGGAACCTCATCAGGAGATGC	25	FAM	BHQ-1
	4	Probe RdRP_SARSr-P1	CCAGGTGGWACRTCATCMGGTGATGC	26	FAM	BHQ-1
E gene	5	Primer E_Sarbeco_F1	ACAGGTACGTTAATAGTTAATAGCGT	26		
	6	Primer E_Sarbeco_R2	ATATTGCAGCAGTACGCACACA	22		
	7	Probe E_Sarbeco_P1	ACACTAGCCATCCTTACTGCGCTTCG	26	FAM	BHQ-1
N gene	8	N_Sarbeco_F1	CACATTGGCACCCGCAATC	19		
	9	N_Sarbeco_R1	GAGGAACGAGAAGAGGCTTG	20		
	10	N_Sarbeco_P1	ACTTCCTCAAGGAACAACATTGCCA	25	FAM	BHQ-1

► Primer/Probe Set 02

Name	Description		Oligonucleotide Sequence (5'>3')	Length	5'mod	3'mod
Target 1 (ORF1ab)	11	Forward Primer 1ab	CCCTGTGGGTTTTACTTAA	21		
	12	Reverse Primer 1ab	ACGATTGTGCATCAGCTGA	19		
	13	Probe 1ab	CCGTCTGCGGTATGTGGAAAGGTTATGG	28	FAM	BHQ-1
Target 2 (N)	14	Forward Primer N	GGGGAAGTTCTCTGCTAGAAT	22		
	15	Reverse Primer N	CAGACATTTTGCTCTCAAGCTG	22		
	16	Probe N	TTGCTGCTGCTTGACAGATT	20	FAM	TAMRA

► Primer/Probe Set 03

Name	Description		Oligonucleotide Sequence (5'>3')	Length	5'mod	3'mod
RP-F	17	RNAse P Forward Primer	AGATTTGACCTGCGAGCG	19		
RP-R	18	RNAse P Reverse Primer	GAGCGGCTGTCTCCACAAGT	20		
RP-P	19	RNAse P Probe	TTCTGACCTGAAGGCTCTGCGCG	23	FAM	BHQ-1
2019-nCoV_N1-F	20	2019-nCoV_N1 Forward Primer	GACCCCAAATCAGCGAAAT	20		
2019-nCoV_N1-R	21	2019-nCoV_N1 Reverse Primer	TCTGGTTACTGCCAGTTGAATCTG	24		
2019-nCoV_N1-P	22	2019-nCoV_N1 Probe	ACCCCGCATTACGTTTGGTGGACC	24	FAM	BHQ-1
2019-nCoV_N2-F	23	2019-nCoV_N2 Forward Primer	TTACAAACATTGGCCGCAAA	20		
2019-nCoV_N2-R	24	2019-nCoV_N2 Reverse Primer	GCGCGACATTCCGAAGAA	18		
2019-nCoV_N2-P	25	2019-nCoV_N2 Probe	ACAATTTGCCCCAGCGCTTCAG	23	FAM	BHQ-1
2019-nCoV_N3-F	26	2019-nCoV_N3 Forward Primer	GGGAGCCTTGAATACACCAAAA	22		
2019-nCoV_N3-R	27	2019-nCoV_N3 Reverse Primer	TGTAGCACGATTGCAGCATTG	21		
2019-nCoV_N3-P	28	2019-nCoV_N3 Probe	AYCACATTGGCACCCGCAATCCTG	24	FAM	BHQ-1

► Set 04 Control

Synbio Technologies has obtained the following viral RNA through *in vitro* transcription and prepared positive controls, which helps to control reverse transcription efficiency, achieve standard unification, reduce false negative rates, make test results more reliable, and greatly improve the accuracy of kit detection.


Name	Description		Transcription Sequence
RdRP gen	29	RdRP standards	GTGAAATGGTCATGTGTGGCGGTTCACTATATGTTAAACCAGG TGGAACCTCATCA GGAGATGCCACAACCTGCTTATGCTAATAGTGTTTTAAACATTTG
E gene	30	E gene standards	ACAGGTACGTTAATAGTTAATAGCGTACTTCTTTTCTTGCTTT CGTGGTA TTCTTGCTAGTTACACTAGCCATCCTTACTGCGCTTCGATTGT GTGCGTACTGCTGCAATAT

CHAPTER 9

COVID-19 DNA Solutions

Name	Description		Transcription Sequence
Target 1 (ORF1ab)	31	ORF1ab standards	CCCTGTGGGTTTTACACTTAAAAACACAGTCTGTACCGTCTGC GGTATGTGGAAAGGTTAT GGCTGTAGTTGTGATCAACTCCGCGAACCCATGCTTCAGTCA GCTGATGCACAATCGT
Target 2 (N)	32	N standards	GGGGAACCTTCTCCTGCTAGAAATGGCTGGCAATGGCGGTGATG CTGCTCTTGCTTTGCTGC TGCTTGACAGATTGAACCAGCTTGAGAGCAAAATGTCTG

Name	Description		Transcription Sequence
Endogenous control	33	RNase P	AGATTTGGACCTGCGAGCGGGTTCTGACCTGAAGGCTCTGCG CGGACTTGTGGAGACAGCCGCTC



* The above figure shows the sample map of COVID-19 gene RNA extracted by Synbio Technologies (N-gene part, ORF1ab part, E-gene part, RdRP gene part).

Sequences above were collected from WHO, US CDC, and China CDC official website documents. Synbio Technologies’s related products have been added with RNase P, which meets the double detection standard and can greatly reduce the “false negative” phenomenon in the nucleic acid detection process and improve the detection accuracy.

9.2 Gene Synthesis & Protein Expression for COVID-19

Synbio Technologies’s Syno® DNA synthesis platform provides fast gene and fragment synthesis, high-throughput synthesis, etc. The independently developed CI sequence difficulty index system and NG™ Codon Optimization Technology can facilitate the accurate, fast, and high-quality synthesis of any sequence.

Early in the pandemic, we completed the synthesis of all structural and non-structural protein genes or gene fragments of COVID-19. This has facilitated the acquisition of viral recombinant proteins and the study of their protein structure biology. In addition, when the new COVID-19 variant Omicron was discovered, we responded quickly and successfully synthesized Omicron’s wild-type Spike protein gene and codon-optimized Spike protein gene. We are committed to providing strong support for downstream vaccine development and neutralizing antibody research.

Competitive Advantages

Speed Guarantee

All COVID-19 related orders will be given production priority with 24/7 professional online support

Production Capacity

Our patented NG™ Codon Optimization Technology optimizes genes at no extra cost to you, and we also have multiple protein expression systems available alongside customized production

Product Supply

Omicron wild type Spike protein gene and optimized Spike protein gene are available in stock

Omicron Product Specifications

Products	Amount	Availability
Wild Type Spike Protein Gene	10 µg	In Stock
Optimized Spike Protein Gene	10 µg	In Stock




Fig.1 Wild type Omicron Spike protein gene




Fig.2 Optimized Omicron Spike protein gene


CHAPTER 10

DNA Data Storage

DNA storage technology (using DNA as a data carrier) is a future-oriented technology with epoch-making significance. It uses the four bases of DNA ("A", "G", "T", "C") arranged in a specific order to form a code for information. As a technology that uses synthetic DNA as a storage medium, it can store text, images, audio and video, and then read subsequently read the data.

Scientists at Synbio Technologies have mastered the next-generation synthetic biology technology platform that could dramatically reduce the manufacturing cost of DNA synthesis. We have developed DNA Studio™ to enable bidirectional transcoding between "A, T, C, G" and digital information. Then, according to the code, the DNA sequences are rapidly and accurately synthesized on a chip in large quantities. We sincerely invite you to join in the development of DNA digital storage systems!


Competitive Advantages



DNA Data Storage

DNA Data Storage Case Study

Synthesized and stored in DNA via our powerful Syno® HT DNA synthesis platform, the Analects of Confucius was transformed into As, Ts, Cs, and Gs using our invented code. With DNA Studio™, the recovery rate after sequencing decoding was 99.9%. This is an example of hard storage, which we published in 2017.



Our R&D department transformed a Chinese poem "The Analects of Confucius" into yeast.

The process is: store-transfer-read-decode. This is a completely soft storage. This poem was copied 10 billion times via this fast and economic method to become the most printed work in the world.

10.1 Your IP, Our DNA

Synbio Technologies has launched the world's first DNA Solution to protect our customers' proprietary information by encrypting and synthesizing it within a sequence of DNA. To accomplish this, Synbio Technologies utilizes our patented DNA Studio™ which accurately synthesizes a tailor-made sequence that acts as an exclusive DNA "security code". By utilizing our Syno® Synthesis Platforms, this unique sequence, or "security code", can be synthesized and placed into any position of the DNA. This allows the security information to be synthesized and stored as a biological "watermark".

How it Works



1. Based on our developed DNA Studio™, Synbio Technologies can provide the requested DNA sequence at any position in the vector as a "security code".
2. With the encryption mechanisms and the Syno® Synthesis Platforms, cryptographic sequences can be synthesized into any section of the requested DNA sequence.
3. The "security code" can also be constructed to any vector requested by the customer.
4. The interpretation process can be easily completed by sequencing and analysis.

DNA Studio™ IP Protection Solution ALLOWS YOU TO



Provide customized **IP protection** for your embedded "security code".

Make your DNA construction and synthesis **"smart"**.